

Multifunctional UV Cabinet

C-65 UV Cabinet

C-65 Chromato-Vue Cabinet accommodates one or two EL Series 8-watt lamps ordered separately. Several lamps are listed below. For additional lamps, see pages 19-20. Use of two lamps illuminates the interior uniformly from two sides of the cabinet.

- The bottom panel of the cabinet is removable
- Includes an interior overhead white light
- Access curtain blocks ambient light
- Soft rubber viewport provides a comfortable viewing area with space for eyeglasses
- Lightweight access curtain blocks out external light from the cabinet interior

C-65 UV Cabinet is designed to accommodate one or two of the EL Series 8-watt Lamps, they can be used with the C-65 cabinet or removed for handheld use (See pages 19-20).



External Dim: 19.75L x 16.5W x 15.75H in. (602 x 419 x 400 mm)
Internal Dim: 17L x 14W x 6H in. (533 x 330 x 400 mm)
Weight: 13.5 lbs.

Order Number

Description

115 V	230 V	Multifunctional UV Cabinet
95-0257-01	95-0257-02	C-65 Mid-Size UV Viewing Cabinet



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.