


Ultraviolet Lamps

Analytik Jena manufactures a large selection of Mineralight, Blak-Ray and UV lamps to accommodate many applications and budgets.

Ultraviolet Lamps: Hand & Specialty

Compact UV Lamps

Compact UV Lamps (4-watts) feature models with longwave (365 nm), shortwave (254 nm), or multiple UV. J-124 Stand for use with Compact UV Lamps.



Dimensions: 7.8L x 2.8W x 2.1H in. (198 x 71 x 53 mm), 1 lb.
J-124 Stand: 6.5H x 7.5W x 4.8D in. (165 x 190 x 122 mm), 1 lb



- Combination wavelength version comes with a snap-on selector plate for switching between wavelengths
- Lightweight, easy-to-handle lamp design
- Comes in 12V style (for cigarette lighter operation)

Order Number		Description
115 V	230 V	Compact UV Lamps, Replacement Tubes and Accessories
95-0016-14	95-0016-15	UVG-11, 4-watt, 254 nm UV
95-0017-09	95-0017-10	UVGL-15 (split-tube), 4-watt, 365 nm/254 nm UV
95-0018-02	95-0018-03	UVL-21, 4-watt, 365 nm UV
95-0021-12	95-0021-10	UVGL-25, 4-watt, 365 nm/254 nm UV
95-0019-01	95-0019-03	UVL-23, 4-watt, 365 nm BLB UV
95-0019-14 (12V)		UVL-23R, 4-watt, 365 nm BL UV
34-0003-01		Tube, 4-watt, 254 nm UV
34-0004-01		Tube, 4-watt, 365/254 nm UV
34-0005-01		Tube, 4-watt, 365 nm UV
34-0010-10		Tube, 4-watt, 365 nm, BLB UV
98-0020-03		J-124 Lamp Stand



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.