



#### Cabinet for chemicals

##### CS-CLASSIC-G model CS.195.105.WDFW:0003

Body colour Anthracite grey (RAL 7016), with wing doors and glass cutout in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30529



#### Function / construction:

- **robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **no unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)
- **easy alignment:** adjusting aids to compensate for uneven floor
- **ventilation:** natural ventilation due to ventilation openings at the cabinet bottom, ready for connection to a technical exhaust system

#### Available equipment

- 2 shelves (height 1105 mm) or 3 shelves (height 1950 mm) galvanised with bottom collecting sump
- doors with glass cutout – the hazardous materials can be seen at a glance, even when the doors are closed

top  
price  
**AVAILABLE  
FROM STOCK  
IMMEDIATELY!**



#### Cabinet for chemicals

##### CS-CLASSIC-G model CS.110.105.WDFW:0002

Body colour Anthracite grey (RAL 7016), with wing doors and glass cutout in Traffic red (RAL 3020), interior equipment with 2 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30527



#### Function / construction:

- **robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **no unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)
- **easy alignment:** adjusting aids to compensate for uneven floor
- **ventilation:** natural ventilation due to ventilation openings at the cabinet bottom, ready for connection to a technical exhaust system

#### Available equipment

- 2 shelves (height 1105 mm) or 3 shelves (height 1950 mm) galvanised with bottom collecting sump

#### Cabinet for chemicals

##### CS-CLASSIC model CS.195.105:0003

Body colour Anthracite grey (RAL 7016), with wing doors in Traffic red (RAL 3020), interior equipment with 3 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30528



top  
price

**AVAILABLE  
FROM STOCK  
IMMEDIATELY!**



#### Cabinet for chemicals

##### CS-CLASSIC model CS.110.105:0002

Body colour Anthracite grey (RAL 7016), with wing doors in Traffic red (RAL 3020), interior equipment with 2 x shelf, 1 x bottom collecting sump (galvanised sheet steel)

Order No. 30526



**Safe and proper storage of non flammable chemicals and toxic substances in work spaces**

**Function / construction:**

- **robust construction and longevity:** robust, high-quality, scratch-resistant outer body with powder-coated textured surface
- **no unauthorised use:** doors lockable with cylinder locking (integrated in the turning handle)
- **easy alignment:** adjusting aids to compensate for uneven floor

- **ventilation:** natural ventilation due to ventilation openings at the cabinet bottom, ready for connection to a technical exhaust system

**Available equipment**

- shelves – optionally in galvanised or painted version, height adjustable
- perforated insert useable as storage area
- bottom collecting sump – optionally in galvanised or painted version, liquid-tight
- pull-out shelves with removable plastic sump

**Cabinet for chemicals**

**CS-CLASSIC model CS.195.105**

Body colour Anthracite grey (RAL 7016), with wing doors in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)

**Order No. 30495-041-30508**



**Cabinet for chemicals**

**CS-CLASSIC-G model CS.195.105.WDFW**

Body colour Anthracite grey (RAL 7016), with wing doors incl. glass cutouts in Pure white (RAL 9010), interior equipment with 4 x pull-out shelf (sheet steel/polyethylene powder-coated/-)

**Order No. 30496-050-30516**



**Cabinet for chemicals**

**CS-PHOENIX model CS.195.105.FD**

Body colour Anthracite grey (RAL 7016), with folding doors in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)

**Order No. 30497-041-30508**



**Cabinet for chemicals**


**CS-CLASSIC model CS.195.105**

Body colour Anthracite grey (RAL 7016), with wing doors in Gentian blue (RAL 5010), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel galvanised)


**Order No. 30495-047-30508**




Front view



Side view



Top view



## Technical data

External dimensions W x D x H  
Internal dimensions W x D x H  
Weight without interior equip.  
Extraction air

## All models

1055 x 520 x 1950 mm  
996 x 493 x 1835 mm  
75 kg  
75 DN

## Technical data

Air change **10 times**  
Diff. pressure (air change **10 times**)  
Maximum load  
Distributed load

## All models

9 m<sup>3</sup>/h  
2 Pa  
600 kg  
421.00 kg/m<sup>2</sup>

7 door colours are available **at no extra cost**:



Order example - how to order correctly



| Model           | Information / equipment                                         | Order No. |
|-----------------|-----------------------------------------------------------------|-----------|
| CS.195.105      | without interior equipment                                      | 30495     |
| CS.195.105.WDFW | without interior equipment, with wing doors incl. glass cutouts | 30496     |
| CS.195.105.FD   | without interior equipment, with folding doors                  | 30497     |

| Body colour              | Door colour             | Order No. |
|--------------------------|-------------------------|-----------|
| Anthracite grey RAL 7016 | Warning yellow RAL 1004 | 040       |
| Anthracite grey RAL 7016 | Light grey RAL 7035     | 041       |
| Anthracite grey RAL 7016 | Traffic red RAL 3020    | 046       |
| Anthracite grey RAL 7016 | Gentian blue RAL 5010   | 047       |
| Anthracite grey RAL 7016 | Reseda green RAL 6011   | 048       |
| Anthracite grey RAL 7016 | asescos silver          | 049       |
| Anthracite grey RAL 7016 | Pure white RAL 9010     | 050       |

| Interior equipment packages (can only be ordered in combination with the cabinet) | Material                                            | Order No. |
|-----------------------------------------------------------------------------------|-----------------------------------------------------|-----------|
| 3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)               | sheet steel galvanised                              | 30508     |
| 3x shelf, 1x bottom collecting sump (V=33.0L)                                     | sheet steel galvanised                              | 30509     |
| 3x shelf, 1x perforated insert, 1x bottom collecting sump (V=33.0L)               | sheet steel powder-coated, sheet steel galvanised   | 30510     |
| 3x shelf, 1x bottom collecting sump (V=33.0L)                                     | sheet steel powder-coated, sheet steel galvanised   | 30511     |
| 3x shelf, 1x bottom collecting sump with grid (V=33.0L)                           | sheet steel galvanised, sheet steel painted blue    | 30512     |
| 3x shelf, 1x bottom collecting sump with grid (V=33.0L)                           | sheet steel powder-coated, sheet steel painted blue | 30513     |
| 3x shelf, 1x bottom collecting sump (large) with grid (V=66.0L)                   | sheet steel galvanised, sheet steel painted blue    | 30514     |
| 3x shelf, 1x bottom collecting sump (large) with grid (V=66.0L)                   | sheet steel powder-coated, sheet steel painted blue | 30515     |
| 4x pull-out shelf (V=33.0L)                                                       | sheet steel/polyethylene powder-coated/- grey/black | 30516     |
| 5x pull-out shelf (V=33.0L)                                                       | sheet steel/polyethylene powder-coated/- grey/black | 30517     |
| 6x pull-out shelf (V=33.0L)                                                       | sheet steel/polyethylene powder-coated/- grey/black | 30518     |

| Interior equipment                       | Material                             | W x D x H (mm)   | 1)     | 2)    | 3)    | Order No. |
|------------------------------------------|--------------------------------------|------------------|--------|-------|-------|-----------|
| shelf                                    | sheet steel galvanised               | 993 x 492 x 30   | 100.00 |       |       | 30519     |
| shelf                                    | sheet steel powder-coated            | 993 x 492 x 30   | 100.00 |       |       | 30520     |
| perforated insert                        | sheet steel galvanised               | 976 x 427 x 41   | 60.00  |       |       | 30521     |
| bottom collecting sump                   | sheet steel galvanised               | 981 x 430 x 85*  |        | 33.00 | 30.00 | 30522     |
| bottom collecting sump with grid         | sheet steel coated blue              | 984 x 424 x 106* |        | 33.00 | 30.00 | 30523     |
| bottom collecting sump (large) with grid | sheet steel coated blue              | 984 x 382 x 175* |        | 66.00 | 60.00 | 30524     |
| pull-out shelf, with PP sump             | sheet steel/polypropylene grey/black | 898 x 418 x 80*  | 30.00  | 34.00 | 30.91 | 30525     |
| box for chemicals                        | PP blue                              | 410 x 300 x 260  |        | 18.00 | 16.36 | 30038     |

| Accessories                                    | Order No. |
|------------------------------------------------|-----------|
| extraction unit without exhaust air monitoring | 14218     |

\* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.