



Safe and approved storage of aggressive and flammable hazardous substances in working areas

Function / construction:

- **robust construction and longevity:** triple hinged door, body made of powder-coated sheet steel, safety elements assembled outside the cabinet for increased protection against corrosion
- **two in one:** 2 compartments (vertical separation) allow the storage of flammable liquids and aggressive hazardous materials in one cabinet
- **user-friendly and safe (folding door compartment):** door open arrest system (standard), convenient and safe removal or storage of containers

- **no unauthorised use:** doors lockable with profile cylinder (integration in an existing locking system possible) / cylinder lock

Available equipment

- Compartment for the storage of flammable liquids
- shelves with load capacity of 75 kg
- tray shelves with load capacity of 150 kg
- perforated insert useable as storage area
- storage box

- Compartment for the storage of aggressive, non flammable hazardous materials
- metal-free tray shelves and collecting sumps

Combi safety storage cabinet

K-PHOENIX-30 model K30.197.120.MV.FWAS

in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated), 4 x tray shelf (PP)

Order No. 30342-001-30344 4

Left compartment



Right compartment



Easy transport

With integrated transport base for an easy in-house transport



Bottom collecting sump

High storage capacity, optimal for the storage of larger containers, safe collection of leakages



Extraction unit

Safe ventilation of safety storage cabinets

plug-in – fast and easy to install

Further information and products from page 397

Order No. 14218 2



Shelf

Height adjustable (32 mm increments), with 30 mm upstand



Drawers

High storage capacity
Ideal for the storage of small containers

**Combi safety storage cabinet****K-PHOENIX-30 model K30.197.120.MV.FWAS**

in Light grey (RAL 7035), interior equipment with 3 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated), 6 x tray shelf (PP)

Order No. 30342-001-30345

**Hazardous substances center****K-PHOENIX-30 model K30.197.120.MC.FWAS**

in Light grey (RAL 7035), interior equipment with 1 x storage box, 2 x shelf, 1 x perforated insert, 1 x bottom collecting sump (sheet steel powder-coated), 4 x tray shelf (PP)

Order No. 30343-001-30758

Order example - how to order correctly

30342



001



30761

**Model**

K30.197.120.MV.FWAS
K30.197.120.MC.FWAS

Information / equipment

without interior equipment, with vertical dividing wall
without interior equipment, hazardous substances center (including lockable storage box on shelf)

Order No.

30342
30343

**Body colour**

Light grey RAL 7035

Door colour

Light grey RAL 7035

Order No.

001

**Interior equipment packages (can only be ordered in combination with the cabinet)****K30.197.120.MV.FWAS**

left: 3x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 4x tray shelf (V=18.0L)
left: 3x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 6x tray shelf (V=18.0L)
left: 3x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 4x tray shelf (V=11.0L)
left: 3x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 6x tray shelf (V=11.0L)

Material

sheet steel powder-coated RAL 7035, PP grey
sheet steel powder-coated RAL 7035, PP grey
stainless steel, PP grey
stainless steel, PP grey

Order No.

30344
30345
30346
30347

**Interior equipment packages (can only be ordered in combination with the cabinet)****K30.197.120.MC.FWAS**

left: 1x storage box, 2x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 4x tray shelf (V=18.0L)
left: 1x storage box, 2x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 6x tray shelf (V=18.0L)
left: 1x storage box, 2x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 4x tray shelf (V=11.0L)
left: 1x storage box, 2x shelf, 1x perforated insert, 1x bottom collecting sump (V=270L)
right: 6x tray shelf (V=11.0L)

Material

sheet steel powder-coated RAL 7035, PP grey
sheet steel powder-coated RAL 7035, PP grey
sheet steel powder-coated, stainless steel, PP grey
sheet steel powder-coated, stainless steel, PP grey

Order No.

30758
30759
30760
30761



K-PHOENIX-30 ◀ 120 cm ▶

Individual configuration of model **K30.197.120.MV(MC).FWAS**

 Model	Information / equipment	W x D x H (mm)	Order No.
K30.197.120.MV.FWAS	without interior equipment, with vertical dividing wall	1202 x 598 x 1968	30342
K30.197.120.MC.FWAS	without interior equipment, hazardous substances center (including lockable storage box on shelf)	1202 x 599 x 1968	30343

 Body colour	Door colour	Order No.
Light grey RAL 7035	Light grey RAL 7035	001 


 Interior equipment individual parts (right comp.)	Material	W x D x H (mm)	 ¹⁾	 ²⁾	 ³⁾	Order No.
 tray shelf	panel/PP melamine resin-coated grey	407 x 468 x 80*	25.00	18.00	16.36	30341 

 Interior equipment individual parts (left comp.)	Material	W x D x H (mm)	 ¹⁾	 ²⁾	 ³⁾	Order No.
 shelf	sheet steel powder-coated RAL 7035	465 x 497 x 40*	75.00			30311 
	stainless steel	465 x 497 x 40*	75.00			30312 
 tray shelf	panel/PP melamine resin-coated grey	407 x 468 x 80*	25.00	18.00	16.36	30341 
 perforated insert	sheet steel powder-coated RAL 7035	420 x 439 x 60	75.00			5555 
	stainless steel	425 x 439 x 60	75.00			5776 
 bottom collecting sump	sheet steel powder-coated RAL 7035	504 x 478 x 120*		2700	24.55	9445 
	stainless steel	504 x 478 x 120*		2700	24.55	10533 
 storage box, left hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30029 
 storage box, right hinged	sheet steel powder-coated RAL 7035	380 x 420 x 360				30030 
 pipe lead-through	sheet steel powder-coated RAL 7035	182 x 104 x 48				30035 
 box for chemicals	PP blue	410 x 300 x 260		18.00	16.36	30038 
 document holder, DIN A5	polystyrene					30793 

 Accessories	Order No.
extraction unit without exhaust air monitoring	14218 

* usable storage space W x D (mm)

1) load capacity (kg), 2) retention capacity (ltr), 3) maximum container volume in the cabinet (ltr)



Technical data


External dimensions W x D x H	mm	1202 x 598 x 1968
Internal dimensions W x D x H	mm	left: 532 x 519 x 1691 right: 486 x 519 x 1691
Weight without interior equip.	kg	330
Air change 10 times	m ³ /h	9
Diff. pressure (air change 10 times EK5/AK4)	Pa	5
Diff. pressure (air change 10 times)	Pa	5
Maximum load	kg	600
Distributed load	kg/m ²	484

Transport base

Entry width transport base	mm	1071
Entry height transport base	mm	100

K30.197.120.MV.FWAS(MC.FWAS)

Top view




K30.197.120.MV(MC).FWAS standard storage heights for shelves (left compartment)

Shelves/number	9	8	7	6	5	4+1	4	3	2
Distance ④ (mm)							365		
Distance ③ (mm)							285	395	
Distance ② (mm)							315	380	
Distance ① (mm)							315	410	

K30.197.120.MV(MC).FWAS standard storage heights for shelves (right compartment)

Shelves/number	9	8	7	6	5	4+1	4	3	2
Distance ⑥ (mm)				170					
Distance ⑤ (mm)				280					
Distance ④ (mm)				280			375		
Distance ③ (mm)				180			380		
Distance ② (mm)				280			380		
Distance ① (mm)				280			380		





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.