

MC-12™

High Speed Microcentrifuge



- + *Less than 9 inches wide*
- + *Large LCD display with rcf conversion*
- + *Whisper quiet operation*
- + *Stops and opens in just 10 seconds*
- + *12 tube capacity (1.5ml) / 16,100 xg*



High Speed Range



Large LED Controls

Benchmark's MC-12 is the first and only compact microcentrifuge capable of achieving centrifugal speeds up to the critical threshold of 16,000 xg, often held as the standard for many molecular biology/DNA protocols. It's brushless, maintenance free drive system makes it powerful, yet compact enough to allow each work station to be fully equipped with a personal centrifuge.

The MC-12 defines a new class of microcentrifuge: compact, quick, powerful and affordable. All aspects of microcentrifuge performance were taken into consideration and no compromises were made. In addition to its unsurpassed speed, the MC-12 also delivers shortened acceleration/deceleration times, whisper quiet operation and its computer-designed airflow pattern minimizes sample temperature increases during longer runs at high speeds.

All operating parameters are set using a single, conveniently located control knob and the large, bright LCD provides a clear indication of time and rpm or g-force. A separate momentary spin button permits instant spin downs, without setting or disturbing any preset parameters. The aircraft grade aluminum rotor accepts 1.5 to 2.0ml tubes directly, as well as 0.5 and 0.2ml tubes with optional adapters. For applications requiring sub-ambient temperatures, the MC-12 is cold room safe.

Technical Data:

Speed:	15,500rpm / 16,100 xg
Capacity:	12 x 1.5/2.0ml
Deceleration:	10 seconds
Acceleration:	12 seconds
Increment:	100 rpm / 100 rcf
Timer:	20 sec. to 99 min. / cont.
Momentary:	Yes
Dimensions:	10.6 x 8.7 x 5.9 in. 27 x 22 x 15 cm
Weight:	10.5 lbs / 4.8 kg
Electrical:	115 to 230V, 50-60Hz, 100W

Ordering Information:

C1612*	MC-12 High Speed Microcentrifuge with 12 x 1.5/2.0ml rotor, 115V
C1008-A5-6	Adapters, 0.5ml, pack of 6
C1008-A2-6	Adapters, 0.2ml, pack of 6

*115V with US plug. To order in 230V, add (-E) to the item number.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.