

SB-12L Shaking Water Bath



- + *Orbital shaking motion*
- + *Included spring platform accepts tubes, bottles and flasks*
- + *Compact design, conserves bench space*
- + *Digital control and display of time, temperature, and speed*

This compact shaking water bath utilizes the well proven, belt-less motor drive system from our popular Orbi-Shaker™ Series, with the addition of precision temperature control up to 80°C.

The seamless, stainless steel chamber has a 12 liter capacity and includes a universal spring-loaded platform for securing 15 & 50ml test tubes, as well as Erlenmeyer flasks and media bottles up to 1 liter. A hinged and removable lid is also included to reduce energy waste, prevent splashing and divert condensation away from samples.

Easy to use controls allow the user to digitally select and monitor time, temperature and shaking speed, while a constant monitoring system ensures accuracy over each parameter.

Optional test tube racks are available for microtubes and custom-designed racks can also be made available by request.

Technical Data:

Temperature range:	Ambient +5°C to 80°C
Capacity:	12 Liters
Temperature accuracy:	+/- 0.3°C at 37°C
Temperature uniformity:	+/- 0.2°C at 37°C
Temperature increment:	0.1°C
Shaking speed range:	30 to 200rpm (increments of 1rpm)
Shaking orbit:	25mm (1")
Timer:	1 min. to 48 hours or continuous
Dimensions:	12.2 x 13.5 x 12 in. 31.7 x 34.3 x 30 cm
Electrical:	115V, 60Hz or 230V, 50Hz / 300W
Weight:	31 lbs / 14kg

Capacity Chart

	Bottles/Tubes	Flasks
0.5mL	40	N/A
1.5/2.0mL	40	N/A
15mL	63	N/A
50mL	15/30	12
125mL	8	8
250mL	4	3
500mL*	2	2
1000mL*	2	1

*Height of vessel may limit ability to close included lid.

Ordering Information

SB0012*	SB-12L Shaking Water Bath with universal spring platform for tubes, bottles and flasks
SB0012-T1520	Optional rack for 1.5/2.0mL micro tubes
SB0012-T5	Optional rack for 0.5mL micro tubes

*115V with US plug, for 230V (EU plug), please add (-E).



SB-12L with 15mL conical tubes



SB-12L with 50mL conical tubes



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.