



Portable Containment Solutions

## PCR Still Air Bubble



The PCR Still Air Bubble provides a comfortable to use, easy to clean enclosed environment for the preparation of PCR reagents. The integrated UV-C Germicidal Lamp provides effective bio-decontamination.

- UV-C Lamp as Standard
- Easy To Clean and Maintain
- Ergonomic Design
- Easy to Install
- Lightweight and Portable
- Pre Tested and Validated



Lab Bubble is a registered European Design: No. 001324610-0001

## FUNCTIONAL SPECIFICATION



The PCR Still Air Bubble is designed to provide an enclosed sample handling area to help prevent contaminants affecting PCR reagent preparation.

Each enclosure is fitted with a UV-C germicidal lamp tested to ensure a proven 254nm wavelength is available to provide bio-decontamination of work surfaces before and after use.

## COLOURS AVAILABLE



## DIMENSIONS



## TECHNICAL SPECIFICATIONS

**Performance:** UV-C Germicidal Lamp (254nm)  
30 minute Exposure

**Power:** 90-240V 1ph 50/60hz Load 1.0A

## CONSTRUCTION

**Weight:** 10kg

**Base:** 10mm Chemically Resistant,  
Process Modified Polypropylene

**Top:** 6mm Chemically Resistant, Impact  
Modified Acrylic

## INSTALLATION

User Install, No Engineer Required  
Pre Validated and Tested System  
Documentation Pack Provided  
No Facility Air Supply Required





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.