

# EXPellPLUS

INDIVIDUALLY WRAPPED TIPS

EXPellPlus individually wrapped tips are ideal for sensitive applications requiring the strictest aseptic conditions.

Their unique individual packing ensures sterility for every single pipetting cycle, maximizing the exceptional performance of EXPellPlus sterile tips. Assuring the greatest purity standards, EXPellPlus individually wrapped tips are free from Endotoxin, DNA, DNase, RNase and PCR inhibitor.

Available in four volume sizes: 10µL standard and extended length, 200µL and 1000-1250µL, EXPellPlus individually wrapped tips are universally compatible with the majority of pipette brands available on the market. Thanks to the unique composition of plastic polymer, the inner surface of EXPellPlus individually wrapped tips is more hydrophobic, which significantly reduces sample adhesion. These low retention properties minimize sample binding inside the tip, making EXPellPlus tips ideal for applications requiring high accuracy and reproducibility.

EXPellPlus individually wrapped tips are sterile. Batch number and expiration date are printed on each individual wrap.



## Ordering Information

Cat. No.	Description
5030012C	ExpellPlus 10µL, pre-sterile, clear, bag, 5x400 pcs. individually wrapped tips
5030042C	ExpellPlus 10µL XL, pre-sterile, long, clear, bag, 5x400 pcs. individually wrapped tips
5030072C	ExpellPlus 200µL, pre-sterile, clear, bag, 5x400 pcs. individually wrapped tips
5130132C	Expell 1000µL (1250µL), pre-sterile, clear, bag, 5x200 pcs. individually wrapped tips





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.