



Display indicative only

## C2004 Small Dual Centrifuge

C2004. (230V 50/60Hz). 1.C2004. (110V 60Hz)

|             |                                 |
|-------------|---------------------------------|
| Speed max   | 500-4,000 Rpm (10 Rpm steps)    |
| Rcf max     | 2,200 G                         |
| Timer       | 0-99 Mins & Hold (30 sec steps) |
| Dims HWD    | 275 x 325 x 470mm               |
| Weight      | 21.8Kg (without rotor)          |
| Power       | 160 watts                       |
| Memory      | 10 programs                     |
| Accel Rates | 10 programs                     |
| Decel Rates | 10 programs                     |

## Swing Out Rotor

Adaptors required see page 17

## Rotors 4000 Rpm



(Requires 4 buckets B5010)

| Rotor                    | BRK5324    | BRK5308    | BRK1000    |
|--------------------------|------------|------------|------------|
| Rotor type               | 24 x 15ml  | 8 x 50ml   | 4 x 100ml  |
| Tube size max            | 17 x 120mm | 30 x 120mm | 45 x 100mm |
| Minimum speed Rpm        | 500        | 500        | 500        |
| Maximum speed Rpm        | 4,000      | 4,000      | 4,000      |
| Maximum rcf (G)          | 2,200      | 2,200      | 2,200      |
| Radius max cms           | 12.3       | 12.3       | 12.3       |
| Sample tube angle °      | 30         | 30         | 0          |
| Acceleration time (secs) | 35         | 35         | 25         |
| Deceleration time (secs) | 35         | 35         | 25         |
| Autoclavable (frequency) | 121°C (20) | 121°C (20) | 121°C (20) |

## Reducers (Pack of 4)



BRK5324



BRK5308



| Rotor         | BRK5324     | BRK5308     |
|---------------|-------------|-------------|
| Part number   | RM05 (5ml)  | RM15 (15ml) |
| Tube size max | 13 x 80mm   | 17 x 120mm  |
| Part number   | RM10 (10ml) | RM25 (25ml) |
| Tube size max | 13 x 100mm  | 25 x 100mm  |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.