

Ceti Steddy Trinocular Stereo Microscope (Diamond-Base)

Product Code: 7500.0000M

Brand: Ceti

Ceti's Steddy Trinocular Stereo Microscope (Diamond-Base) is one of Ceti's most popular zoom microscopes. It provides high quality, wide field of view images, a wide zoom range and is highly robust and extremely well built, suitable for the most demanding uses. This microscope features a convenient 0.65x - 4.5x zoom ratio and a pair of 10X super wide-field (FN22) eyepieces, giving a total magnification range of 6.5 to 45x. The optical quality is extremely good giving a bright, crisp 3D image of the specimen.

It also has dual LED illumination and a photo port to allow the connection of a digital camera to capture images or videos.

This model is one of our most versatile stereo zoom microscopes and is an excellent, professional grade instrument choice for industry, hobbyists, researchers, and educational establishments.



Features

- Diamond base
- Super wide-field 10X, FN22 eyepieces for maximum field of view
- Wide magnification zoom range of 3.35x to 180x
- A trinocular head with photo port
- Dual LED illumination

Specifications

Optical Head	Trinocular, inclined at 45° with photo port
	Interpupillary Distance Adjustment (Siedentopf Type): 55 to 75mm
	Diotropic adjustment on both eyepiece tubes
Eyepieces	10x/22mm super wide-field
Magnification	Zoom Ratio: 0.67 - 45x
	Standard Magnification Range: 6.7 - 45x
	Full Magnification Range (with optional lenses and eyepieces): 3.35 - 180x
Working Distance	165mm
Focusing Mount	Rectangular support with dual coarse focus knobs
Stage	Diamond base with 5W LED illumination
	8W LED ring light
Focusing Knobs	Coaxial coarse with safety stop and tension adjustment
Mains Power	220 - 240V/50 - 60Hz
Supplied With	Two rubber eyecups, dust cover, instruction manual, and blue filter



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.