

Ceti Magnum Trinocular Compound Microscope (10x, 20x, 40x, 60x Polarisation)

Product Code: 2753.0000M

Brand: Ceti

Ceti's Magnum Trinocular Compound Microscope (4x, 10x, 40x, 60x Polarisation) is a popular, economically priced, biological upright microscope. Providing excellent quality images via the wide-field FN20mm eyepieces, this microscope includes four Infinity Corrected Plan Achromatic objectives in the reversed quintuple nosepiece.

It features a robust, T-shaped base, designed to offer maximum stability. Modern, ergonomic features such as a 30° inclined eyepiece tube, ensure a comfortable viewing position over a prolonged period.

This model can also be used in Brightfield mode so it is extremely flexible. Bright, white LED transmitted lighting, with intensity control, is also featured.



Features

- Sturdy base for comfortable viewing over long periods
- Wide-field 10X, FN20 eyepieces with dioptic adjustment
- Quintuple, reversed nosepiece with four Infinity Corrected Plan Achromatic objectives
- Trinocular eyepiece tube as standard to accept digital microscope cameras
- Bright, white LED illumination with intensity control
- Full Köhler illumination

Specifications

Optical Head	Trinocular with photo port, inclined at 30°
	Interpupillary Distance Adjustment (Siedentopf Type): 55 to 75mm
	Diotropic adjustment on one eyepiece
	Analyser
Objectives	4X, 10x, 40x, 60x Infinity Corrected Plan Achromatic
Condenser	Pre-condenser with field diaphragm and polariser
	Focusable and centrable, iris diaphragm
Stage	Round, centrable scale with two stage clips
Focusing Knobs	Coaxial coarse and fine with stop
	Tension adjustment ring
Mains Power	220 - 240V/50 - 60Hz
Illumination, Transmitted	3W LED with intensity control
Supplied With	Blue and green filter, 2 spare fuses, immersion oil (10ml), soft rubber eyecups, dust cover, and instruction manual



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.