

Ceti Varizoom 10 Binocular Stereo Microscope

Product Code: 7700.0000M

Brand: Ceti

The Ceti Varizoom 10 Binocular Stereo Microscope is one of our highest specification stereoscopic zoom microscopes with a parallel zoom system utilising optics of the highest quality.

A sturdy, rectangular focusing mount enables stability whilst zooming and focusing and the broad stand provides a comfortable base when manipulating samples.

Integrated, directional episcopic LED lighting brightly illuminates the sample without the need for awkward goose neck light guides and the transmitted LED lighting provides cool lighting, especially important when working with biological or heat-sensitive samples. Separate intensity controls are located at the rear of the stand so they do not encroach into the working area.

Features

- Broad, flat base for comfortable working
- Super wide-field 10X, FN22 eyepieces for maximum field of view
- Wide magnification zoom range of 4x to 320x utilising parallel zoom optics
- Optional eyepieces and objectives to easily adjust the zoom magnification range
- Dual LED illumination with separate intensity controls



Specifications

Optical Head	Binocular, inclined at 20°
	Interpupillary Distance Adjustment (Sidentopf type): 55 to 75mm
	Dioptric adjustment on both eyepiece tubes
Eyepieces	10x/20mm super wide-field
Zoom Ratio	10:1
Magnification Range (standard)	8x to 80x
Magnification Range (total)	4x to 320x
Working Distance	78mm (with Plan Apo 1x/WF)
Focusing Mount	Rectangular support with dual coarse and fine focusing knobs
Stage	Frosted glass/black and white contrast plate
	2 slide clip holders
Illumination	Directional epi-illumination, attached under the optical head, LED 5W
	Fixed transmitted illumination, built into the base, LED 5W
Mains Power	220 - 240V/50 - 60 Hz
Supplied With	Two rubber eyecups, dust cover, instruction manual, and one blue filter



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.