




omniPAGE Mini Modular Systems

The omniPAGE range of Modular Vertical Gel Systems allow multiple electrophoresis techniques to be performed in the same unit. Using the same main tank and lid, three different modules are interchangeable for PAGE, tube gel and electroblotting techniques.

Each Modular System comprises the Gel Chamber, PAGE Internal Running Module, Casting Base and various accessories. Through the addition of further modules, the packages flexibility is increased to allow Western blotting and Capillary IEF.




1

Capillary IEF module, VS10DCI




2

Western Blotting module, VS10BI



3

Hi-Intensity Blotting module, VS10BI-HI




4

2x Blot Module available, see ordering information below

ACCESSORIES



glass plates with bonded spacers



combs



cool packs

ORDERING INFORMATION


| | | | |
|-------------|---|-----------|--|
| CVS10CES | Complete Mini (10x10cm) Vertical Electrophoresis Modular System , comprising: 1x Mini Vertical Unit, CVS10DSYS (1) which includes: PAGE Module, 2x2mm thick notched glass plates, 2x2mm thick plain glass plates with 1mm thick bonded spacers, 1x blanking plate, 2x combs (1mm thick, 12 samples), 1x casting base, silicone mat, cooling pack plus: 1x Capillary Electrophoresis Module, VS10DCI (2) and 1x Electroblotting Module, VS10BI (3) comprising: internal electroblotting module, 4x compression cassettes for gel sizes up to 10x10cm and 8x fibre pads | | |
| CVS10C2DS | Complete Mini (10 x 10cm) Vertical Electrophoresis & 2-D System , comprising: 1x Mini Vertical Unit, CVS10DSYS (1) which includes: PAGE Module, 2x2mm thick notched glass plates, 2x2mm thick plain glass plates with 1mm thick bonded spacers, 1x blanking plate, 2x combs (1mm thick, 12 samples), 1x casting base, silicone mat, cooling pack plus: 1x Capillary Electrophoresis Module, VS10DCI (2) | | |
| CVS10CBS | Complete Mini (10 x 10cm) Vertical Electrophoresis & Blotting System , comprising: 1x Mini Vertical Unit, CVS10DSYS (1) which includes: PAGE Module, 2x2mm thick notched glass plates, 2x2mm thick plain glass plates with 1mm thick bonded spacers, 1x blanking plate, 2x combs (1mm thick, 12 samples), 1x casting base, silicone mat, cooling pack plus 1x Standard Electroblotting Module, VS10BI (3) | | |
| CVS10CBS-HI | Complete Mini (10 x 10cm) Vertical Electrophoresis & High Intensity Blotting System , comprising: 1x Mini Vertical Unit, CVS10DSYS (1) which includes: PAGE Module, 2x2mm thick notched glass plates, 2x2mm thick plain glass plates with 1mm thick bonded spacers, 1x blanking plate, 2x combs (1mm thick, 12 samples), 1x casting base, silicone mat, cooling pack plus: 1x High Intensity Electroblotting Module, VS10BI-HI (4) | | |
| VS10DCI | omniPAGE Mini Tube Unit (2) | VS10BI-HI | High Intensity omniPAGE Blot Mini insert - includes 2 cassettes and 8 fibre pads (4) |
| VS10BI | OmniBlot Mini Insert - including 4 cassettes, 16 foam pads (3) | SB10 | omniBlot Mini 10 x 10cm Blotting System, 4x cassettes |
| VS10BI-2 | OmniBlot Mini Insert - including 2 cassettes, 8 foam pads | SB10-2 | omniBlot Mini 10 x 10cm Blotting System, 2x cassettes |



OmniPAGE Mini Wide Systems

The Mini Wide vertical gel unit, with a gel width of 20cm, effectively allows double the number of samples to be resolved as the mini unit. This allows consistency of sample comparison on a single gel and is designed for those with greater than 20 samples to compare and resolve. Simple set up using ultra soft silicone seals guarantees trouble free glass plate loading and gel casting.

Each Modular System comprises the Gel Chamber, PAGE Internal Running Module, Casting Base and various accessories. Through the addition of further modules, the packages flexibility is increased to allow Western blotting and Capillary IEF.




Capillary IEF module, VS10WDCI




Western Blotting module, VS10WBI

ACCESSORIES



glass plates with bonded spacers



combs



cool packs

ORDERING INFORMATION

| | | | |
|------------------|--|-------------|--|
| VS10WD | Mini Wide, 20 x 10cm Dual , 2 sets of Glass Plates with 1mm thick bonded Spacers, 2 x 24 sample, 1mm thick combs, cooling pack | | |
| VS10WDSYS | Mini Wide, 20 x 10cm Dual , 2 sets of Glass Plates with 1mm thick bonded Spacers, 2 x 24 sample, 1mm thick combs, cooling pack including caster | | |
| VS20CAST | 20 x 10cm Casting Base | VS10WNGS1 | 20 x 10cm Notched Glass Plates with 1mm Bonded Spacers (pk/2) |
| VS20DCASTM | Replacement Silicone Mat for 20 x 10cm Casting Base | VS10WPGS1 | 20 x 10cm Plain Glass Plates with 1mm Bonded Spacers (pk/2) |
| VS10WDIRM | Inner Running Module | VS10WPGS1.5 | 20 x 10cm Plain Glass Plates with 1.5mm Bonded Spacers (pk/2) |
| VS20-x -LG | Loading guides for RigRunner V-MINI combs, x = comb well number | VS10WPGS2 | 20 x 10cm Plain Glass Plates with 2mm Bonded Spacers (pk/2) |
| VS10WNG | 20 x 10cm Notched Glass Plates 4mm thick (pk/2) | VS10WDP | Blanking Plate, 20 x 10cm |
| VS10WPG | 20 x 10cm Plain Glass Plates 4mm thick (pk/2) | RPW-0.2100 | Replacement Platinum Wire - 0.2mm, 50cm |
| VS10WNGS0.75 | 20 x 10cm Notched Glass Plates with 0.75mm Bonded Spacers (pk/2) | VS20ICB | Cooling Pack |
| VS10WPGS0.75 | 20 x 10cm Plain Glass Plates with 0.75mm Bonded Spacers (pk/2) | | |
| SB10W | Mini Wide Blot Unit, 20 x 10cm System including tank and lid, 4 cassettes, 8 fibre pads, cooling pack | VS10WDC | Mini Wide Tube Gel Unit, 20x10cm with tank and lid, glass capillary tubes, blanking ports and cooling pack |
| VS10WBI | Mini Wide Blot Module - includes 4 cassettes and 8 fibre pads | VS10WDCI | Mini Wide Tube Gel Module - includes glass tubes and blanking ports |
| SB10WC | Mini Wide Blot Cassette | MCT10 | Mini Capillary Tubes, pk/100 |
| SB10WF | Fibre pads - pk/8 | MCT101.5 | Mini Capillary Tubes, 1.5mm, pk/100 |
| VS10WCES | Complete Mini Wide (20x10cm) Vertical Electrophoresis Modular System , comprising: 1x Mini Wide Vertical Unit, VS10WDSYS (1) which includes: PAGE Module, 2x4mm thick notched glass plates, 2x4mm thick plain glass plates with 1mm thick bonded spacers, 1x blanking plate, 2x combs (1mm thick, 24 samples), 1x casting base, silicone mat, cooling pack plus: 1x Wide Electroblothing Module, VS10WBI (3) 1x Wide Capillary Electrophoresis Module, VS10WDCI (2) | | |
| VS10WCBS | Complete Mini Wide (20 x 10cm) Vertical Electrophoresis & Blotting System , comprising: 1x Mini Wide Vertical Unit, VS10WDSYS (1) which includes: 2x4mm thick notched glass plates, 2x4mm thick plain glass plates with 1mm thick bonded spacers, 1x blanking plate, 2x combs (1mm thick, 24 samples), 1x casting base, silicone mat, cooling pack plus: 1x Wide Electroblothing Module, VS10WBI (3) | | |
| VS10WC2DS | Complete Mini Wide (20x10cm) 2-D System , comprising: 1x Mini Wide Vertical Unit, VS10WDSYS (1) which includes: 2x4mm thick notched glass plates, 2x4mm thick plain glass plates with 1mm thick bonded spacers, 1x blanking plate, 2x combs (1mm thick, 24 samples), 1x casting base, silicone mat, cooling pack plus: 1x Wide Capillary Electrophoresis Module, VS10WDCI (3) | | |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.