

NE2D Series Unstirred Digital Water Baths



The **NE2D series** provides a stable temperature environment ideal for meeting quality control requirements, featuring digital PID temperature control for precision. Controls feature auto alarm setting, illuminated on/off switch, heater and over temperature indications.

Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance.

	NE2 - 4D	NE2 - 8D	NE2 - 10D	NE2 - 14D
GENERAL				
Temp. Range	Ambient + 5°C - 100°C			
Temp. Alarm	Over temperature alarm is automatically set 4°C above set point			
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C			
Sensitivity/Uniformity	±0.2°C / ±0.1°C			
Heater Power (220-240V)	400W	800W		1000W
Voltage	120V or 230V models available 50-60Hz			
DIMENSIONS				
Tank Capacity	4 Litre	8 Litre	10 Litre	14 Litre
External (mm)	332w x 185d x 290h	332w x 270d x 290h	332w x 270d x 290h	361w x 332d x 290h
Internal (mm)	300w x 150d x 150h	300w x 240d x 150h	300w x 240d x 200h	325w x 300d x 150h
Water Level Max. (mm)	130		180	130
Water Level Min. (mm)	55			

	NE2 - 22D	NE2 - 28D	NE2 - 56D
GENERAL			
Temp. Range	Ambient + 5°C - 100°C		
Temp. Alarm	Over temperature alarm is automatically set 4°C above set point		
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C		
Sensitivity/Uniformity	±0.2°C / ±0.1°C		
Heater Power (220-240V)	1500W		2500W
Drain (12mm Bore)	Drain Outlet		
Voltage	120V or 230V models available 50-60Hz		
DIMENSIONS			
Tank Capacity	22 Litre	28 Litre	56 Litre
External (mm)	538w x 332d x 290h	538w x 332d x 290h	660w x 540d x 290h
Internal (mm)	500w x 300d x 150h	500w x 300d x 200h	620w x 500d x 200h
Water Level Max. (mm)	130	180	
Water Level Min. (mm)	55		

Accessories on Page 14-17



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.