



NE9-15S

NE9-32S

GENERAL		
Temp. Range	Ambient +5°C to 300°C	
Temp. Alarm	Over temperature alarm is automatically set 5°C above set point	
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C	
Temp. Sensor	pt100 with control diagnostics at power up	
Temp. Fluctuation 150°C	±0.4°C	
Temp. Variation 150°C	±1.8°C	
Temp. Heat up 150°C	25 Minutes	
Temp. Door open recovery 150°C	5 Minutes	
Temp. Control & Run Back Timer	Visual indicator and end messaging 00:01 to 99:59 (hh:mm) activates at set temp	
Temp. Control & Delay Start Timer	Heating commences at a set time in the future, within 24 hours	
Chamber convection	Gravity convection 0% variable to 100% fan speed for forced air circulation	
Exhaust Air Flap	Adjustable 0 to 100% chamber venting	
Access Port	35mm port size, perforated flap cover helps to reduce losses with external probes fitted	
Doors	Single insulated door, lockable with 2 point door engagement for an even seal	
Door Alarm	Door open, heating interrupted, delay before door sounder chimes until door closed, heating resumes	
Shelves fitted & max	1 wire / 3 max	
Shelving Locations	Available every 20mm vertically	
Shelf Max. Load	10kg	20kg
Heater Power	1050W	1100W
Voltage	230V 50/60Hz	
DIMENSIONS		
Chamber Capacity	15 Litre	32 Litre
External Dimensions (mm)	520w x 390d x 600h	600w x 460d x 680h
Internal Dimensions (mm)	320w x 180d x 240h	400w x 250d x 320h
Stackable	Yes	
Feet/Castors	Levelling Feet	

The **NE9 series** Drying Oven, for gravity convection and variable fan speed for forced air circulation and constructed with a stainless steel chamber. Dial controls for Air Flap and Fan Speed settings, bold illuminated LED temperature display offering precise time and temperature control from ambient +5°C to 300°C. The oven can run at set temperature immediately, use run back timer or delay start feature. With energy efficiency in mind, particular attention is paid to providing a good tight door closure and hinge design to maintain a sealed chamber in use. Because the unit has a sealed chamber, less energy is used to hold set temperature and the intelligent controller will pulse energy to maintain set point throughout operation.

NE9-56S		NE9-112S		NE9-240S	
GENERAL					
Temp. Range	Ambient +5°C to 300°C				
Temp. Alarm	Over temperature alarm is automatically set 5°C above set point				
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C				
Temp. Sensor	pt100 with control diagnostics at power up				
Temp. Fluctuation 150°C	±0.4°C				
Temp. Variation 150°C	±1.8°C				
Temp. Heat up 150°C	25 Minutes				
Temp. Door open recovery 150°C	5 Minutes				
Temp. Control & Run Back Timer	Visual indicator and end messaging 00:01 to 99:59 (hh:mm) activates at set temp				
Temp. Control & Delay Start Timer	Heating commences at a set time in the future, within 24 hours				
Chamber convection	Gravity convection 0% variable to 100% fan speed for forced air circulation				
Exhaust Air Flap	Adjustable 0 to 100% chamber venting				
Access Port	35mm port size, perforated flap cover helps to reduce losses with external probes fitted				
Doors	Single insulated door, lockable with 2 point door engagement for an even seal				
Door Alarm	Door open, heating interrupted, delay before door sounder chimes until door closed, heating resumes				
Shelves fitted & max	2 wire / 5 max	2 wire / 8 max		3 wire / 10 max	
Shelving Locations	Available every 20mm vertically				
Shelf Max. Load	20kg				
Heater Power	2150W			3200W	
Voltage	230V 50/60Hz			230V 50/60Hz Blue Power plug 2P+E	
DIMENSIONS					
Chamber Capacity	56 Litre	112 Litre		240 Litre	
External Dimensions (mm)	600w x 605d x 720h	600w x 660d x 980h		800w x 720d x 1160h	
Internal Dimensions (mm)	400w x 400d x 360h	400w x 450d x 620h		600w x 510d x 800h	
Stackable	Yes			No	
Feet/Castors	Levelling Feet				

NE9-400S

NE9-750S

GENERAL		
Temp. Range	Ambient +5°C to 300°C	
Temp. Alarm	Over temperature alarm is automatically set 5°C above set point	
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C	
Temp. Sensor	pt100 with control diagnostics at power up	
Temp. Fluctuation 150°C	±0.8°C	
Temp. Variation 150°C	±2.5°C	
Temp. Heat up 150°C	35 Minutes	
Temp. Door open recovery 150°C	10 Minutes	
Temp. Control & Run Back Timer	Visual indicator and end messaging 00:01 to 99:59 (hh:mm) activates at set temp	
Temp. Control & Delay Start Timer	Heating commences at a set time in the future, within 24 hours	
Chamber convection	Gravity convection 0% variable to 100% fan speed for forced air circulation	
Exhaust Air Flap	Adjustable 0 to 100% chamber venting	
Access Port	35mm port size, perforated flap cover helps to reduce losses with external probes fitted	
Doors	Double door, 2 point locking, even seal	
Door Alarm	Door open, heating interrupted, delay before door sounder chimes until door closed, heating resumes	
Shelves fitted & max	3 wire / 14 max	5 wire / 16 max
Shelving Locations	Available every 20mm vertically	
Shelf Max. Load	20kg	
Heater Power	4250W	5250W
Voltage	230V 50/60Hz Blue Power plug 2P+E	
DIMENSIONS		
Chamber Capacity	400 Litre	750 Litre
External Dimensions (mm)	1000w x 720d x 1470h mm	1240w x 810d x 1625h mm
Internal Dimensions (mm)	800w x 510d x 1040h mm	1040w x 600d x 1200h mm
Stackable	No	
Feet/Castors	Levelling castor	

ACCESSORIES

	Wire Shelves	Perforated Shelves	Additional program function	Factory Calibration 1 point	Temperature Survey 9 point	Temperature Survey 27 point
NE9-15S	WS-15	SP-15	RT-3	CAL	CAL-9	CAL-27
NE9-32S	WS-32	SP-32				
NE9-56S	WS-56	SP-56				
NE9-112S	WS-112	SP-112				
NE9-240S	WS-240	SP-240				
NE9-400S	WS-400	SP-400				
NE9-750S	WS-750	SP-750				



NE9-15S

NE9-32S

GENERAL		
Temp. Range	Ambient +5°C to 300°C	
Temp. Alarm	Over temperature alarm is automatically set 5°C above set point	
Temp. Control	PID Digital, LED display showing actual or set temperature, resolution 0.1°C	
Temp. Sensor	pt100 with control diagnostics at power up	
Temp. Fluctuation 150°C	±0.4°C	
Temp. Variation 150°C	±1.8°C	
Temp. Heat up 150°C	25 Minutes	
Temp. Door open recovery 150°C	5 Minutes	
Temp. Control & Run Back Timer	Visual indicator and end messaging 00:01 to 99:59 (hh:mm) activates at set temp	
Temp. Control & Delay Start Timer	Heating commences at a set time in the future, within 24 hours	
Chamber convection	Gravity convection 0% variable to 100% fan speed for forced air circulation	
Exhaust Air Flap	Adjustable 0 to 100% chamber venting	
Access Port	35mm port size, perforated flap cover helps to reduce losses with external probes fitted	
Doors	Single insulated door, lockable with 2 point door engagement for an even seal	
Door Alarm	Door open, heating interrupted, delay before door sounder chimes until door closed, heating resumes	
Shelves fitted & max	1 wire / 3 max	
Shelving Locations	Available every 20mm vertically	
Shelf Max. Load	10kg	20kg
Heater Power	1050W	1100W
Voltage	230V 50/60Hz	
DIMENSIONS		
Chamber Capacity	15 Litre	32 Litre
External Dimensions (mm)	520w x 390d x 600h	600w x 460d x 680h
Internal Dimensions (mm)	320w x 180d x 240h	400w x 250d x 320h
Stackable	Yes	
Feet/Castors	Levelling Feet	



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.