

NE1D Series Unstirred Digital Water Baths



NE1D-14



Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance. The NE1D series, Clifton Digital Baths are designed for immersion of flasks/bottles/racks in a highly stable temperature environment throughout the water bath.

NE1D-4		NE1D-8		NE1D-14		NE1D-22			
GENERAL									
Temp. Range		Ambient + 5°C - 99°C							
Temp. Control		Digital, LED display showing actual or set temperature, resolution 1°C							
Sensitivity/Uniformity		±0.2°C / ±0.1°C							
Heater Power (220-240V)		400W		800W		1000W		1500W	
Voltage		120V or 230V models available 50-60Hz							
Drain (12mm bore)								Drain Outlet	
DIMENSIONS									
Tank Capacity		4 Litre		8 Litre		14 Litre		22 Litre	
External (mm)		332w x 185d x 290h		332w x 270d x 290h		361w x 332d x 290h		538w x 332d x 290h	
Internal (mm)		300w x 150d x 150h		300w x 240d x 150h		325w x 300d x 150h		500w x 300d x 150h	
Water Level Max. (mm)		130							
Water Level Min. (mm)		55							

Accessories on Page 14-17

Unstirred Bath Accessories - Lids

	Stainless Steel Gabled Lid with Insulated Handle	Stainless Steel Hinged Gable Lid with Insulated Handle	Stainless Steel Flat One Piece Lid with Insulated Handle	Clear Lift Off Gable Lid
UNSTIRRED BATHS				
NE1-2.5	Not Available	Not Available	LD-2.5	Not Available
NE1-2/2/2, NE1D-2/2/2, NE2-2D	SL1-4	Not Available	LD-4	GL1-4AF
NE1-4, NE1D-4, NE1-4/4/4, NE1D- 4/4/4, NE2-4D, NE2D-4/4	SL1-4	Not Available	LD-4	GL1-4AF
NE2-5D	SL1-14	Not Available	LD-14	GL1-14AF
NE1-8, NE1D-8, NE2-8D, NE2D-4/8, NE2D-8/22	SL1-8	Not Available	LD-8	GL1-8AF
NE1-9, NE2-9D	SL1-22	SL1-22H	LD-22	GL1-22AF
NE2-10D	SL1-8	Not Available	LD-8	GL1-8AF
NE1-14, NE1D-14 NE2-14D	SL1-14	SL1-14H	LD-14	GL1-14AF
NE1B-14, NE1B-18W	Not Available	SL1-14H/SS	LD-14	GL1-14AF
NE1-22, NE1D-22, NE2-22D, NE2D-8/22	SL1-22	SL1-22H	LD-22	GL1-22AF
NE1-28, NE2-28D	SL1-22	SL1-22H	LD-22	GL1-22AF
NE1-56, NE2-56D	SL1-56	Not Available	LD-56	Not Available
NE1B-28	SL1-22SS	Not Available	LD-22	GL1-22AF



Stainless Steel Lid with Concentric Rings SLR1-4 : 2 x 105 mm holes SLR1-8 : 4 x 83mm holes SLR2-8 : 2 x 105 mm holes SLR1-14 : 4 x 105mm holes SLR2-14 : 6 x 83mm holes SLR3-14 : 4 x 114.5mm holes SLR1-22 : 6 x 105mm holes SLR3-22 : 6 x 114.5mm holes	Clear Hinged Gable Lid	Polypropylene Spheres: Pack of 200
		
Not Available	Not Available	BP0368 : 1 Pack
SLR1-4	GL1-4	BP0368 : 1 Pack
SLR1-4	GL1-4	BP0368 : 1 Pack
SLR1-14, SLR2-14 and SLR3-14	GL1-14	BP0368 : 2 Packs
SLR1-8 and SLR2-8	GL1-8	BP0368 : 2 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 2 Packs
SLR1-8 and SLR2-8	GL1-8	BP0368 : 2 Packs
SLR1-14, SLR2-14 and SLR3-14	GL1-14	BP0368 : 2 Packs
SLR1-14, SLR2-14 and SLR3-14	GL1-14	BP0368 : 2 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 4 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 4 Packs
Not Available	Not Available	BP0368 : 4 Packs
SLR1-22 and SLR3-22	GL1-22	BP0368 : 4 Packs

Accessories for the NE1B-18W Boiling Bath

ST-1	Raise and Lower Sterilising Tray
SA01409	Pair of Guide Rail Assemblies for ST-1 Tray
SL1-14H/SS	Stainless Steel Hinged Gable Lid to suit 14 litre and 18 litre Boiling Baths
GL1-14AF	Lift Off Gable Lid to suit 14 and 18 litre Boiling Baths
GL1-14	Clear Hinged Gable Lid to suit 14 and 18 litre Boiling Baths

Unstirred Bath Accessories

	Stainless Steel Test Tube Racks	PVC Churchill Racks	Drain Tap
	6870 : 26 holes x 17mm DIA 6871 : 16 holes x 26mm DIA 6872 : 36 holes x 13mm DIA 6873 : 18 holes x 19mm DIA 6900 : 12 holes x 32mm DIA	C18-4 : 28 holes x 18mm DIA C18-8 : 44 holes x 18mm DIA C18-14 : 68 holes x 18mm DIA C25-4 : 28 holes x 25mm DIA C25-8 : 44 holes x 25mm DIA C25-14 : 68 holes x 25mm DIA	
UNSTIRRED BATHS			
NE1-2.5	1 Rack From Above Range	Not Available	Not Available
NE1-2/2/2, NE1D-2/2/2, NE2-2D	1 Rack From Above Range Chamber	Not Available	Not Available
NE1-4, NE1D-4, NE1-4/4/4, NE1D-4/4/4, NE2-4D, NE2D-4/4	1 Rack From Above Range	C18-4 and C25-4	Not Available
NE2-5D	4 Racks From Above Range	Not Available	Not Available
NE1-8, NE1D-8 NE2-8D, NE2D-4/8, NE2D-8/22	2 Racks From Above Range	C18-4 and C25-8	Not Available
NE1-9, NE2-9D	6 Racks From Above Range	Not Available	Not Available
NE2-10D	2 Racks From Above Range	Not Available	Not Available
NE1-14, NE1D-14 NE2-14D	4 Racks From Above Range	C18-14 and C25-14	Not Available
NE1B-14	4 Racks From Above Range	Not Available	Not Available
NE1B-18W	4 Racks From Above Range	Not Available	BX0953
NE1-22, NE1D-22, NE2-22D, NE2D-8/22	6 Racks From Above Range	Not Available	BX0953
NE1B-28, NE1-28, NE2-28D	6 Racks From Above Range	Not Available	BX0953
NE1-56, NE2-56D	12 Racks From Above Range	Not Available	BX0953

Flat n Fold™ Test Tube Racks 6770 : 26 holes x 17mm DIA 6771 : 16 holes x 26mm DIA 6772 : 36 holes x 13mm DIA 6773 : 18 holes x 19mm DIA 6800 : 14 holes x 32mm DIA	Stainless Steel Wire Baskets	Bath Armor™	Thermometer Clip With Bent Stem Thermometer
			
1 Rack From Above Range	Not Available	BA-5 : 4 Litre Pack	TC-1
1 Rack From Above Range Chamber	Not Available	BA-5 : 4 Litre Pack	TC-1
1 Rack From Above Range	BSK-4	BA-5 : 4 Litre Pack	TC-1
4 Racks From Above Range	Not Available	BA-12 : 12 Litre Pack	TC-1
2 Racks From Above Range	BSK-8	BA-12 : 12 Litre Pack	TC-1
6 Racks From Above Range	Not Available	BA-12 : 12 Litre Pack	TC-1
2 Racks From Above Range	Not Available	BA-12 : 12 Litre Pack	TC-1
4 Racks From Above Range	BSK-14	BA-12 : 12 Litre Pack	TC-1
4 Racks From Above Range	BSK-14	Not Available	TC-1
4 Racks From Above Range	Not Available	Not Available	TC-1
6 Racks From Above Range	BSK-22	BA-12 : 2 x 12 Litre Pack Required	TC-1
6 Racks From Above Range	Not Available	BA-12 : 2 x 12 Litre Pack Required	TC-1
12 Racks From Above Range	Not Available	BA-12 : 2 x 12 Litre Pack Required	TC-1

Common Accessories

BX0616	Drain Syphon
BKF-1	Bar Keepers Friend



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.