

Cole-Parmer®

Cole-Parmer® Block Heater, 2 Block, Analogue BH-200

- Analogue control
- Economical price
- 130°C maximum temperature
- Hot warning light to warn that the blocks are too hot to touch
- Uniform and stable temperatures



Cole-Parmer® Block Heater, 2 Block, Analogue BH-200

Designed for the precise heating of test-tubes, microcentrifuge tubes, cuvettes and microtitre plates. With simple analogue setting of temperature for ease of use. “Hot” warning light will flash whenever the block temperature is above 50°C.



Technical Specification

Specification	BH-200
Temperature range	Ambient +8°C to 130°C
Temperature stability at 37°C	±0.1°C
Uniformity within block at 37°C	±0.1°C
Uniformity within block at 130°C	±1°C
Dimensions, mm (w x d x h)	235 x 280 x 115
Net weight, kg (without blocks)	2.1
Heater power	300W
Electrical supply	230V, 50Hz
IP Rating	31

Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Block heater, 2 block, analogue, 130°C	36610-01	BH-200	BH-200-2	SBH130
Block heater analogue Dual 130C 120V	36610-00	BH-200	BH-200-2-120	SBH130/120V/60



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.