

# Cole-Parmer VP-200 Series Vacuum/Pressure Pump – 12 VDC

- Get better results with minimal dead space
- Get long, reliable operation—continuous operation up to 18 psi for 4000 service hours
- Provides pumping with minimal stress and heat buildup

**Note:** Do not use metal fittings.

WARNING: This product can expose you to chemicals including BPA, lead, and phthalates, which is known to the State of California to cause cancer and reproductive harm. For more information, visit [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



Vacuum/Pressure Pump 12 VDC  
07532-25

Free air capacity	Maximum vacuum	Maximum pressure*	Port connections	Motor	Operating conditions	Dimensions	Power	Item number
17.2 L/min (0.61 cfm)	508 mm Hg (20 in Hg)	18 psig (1.24 bar)	3/8" NPT	TEFC	0 to 40 °C (32 to 104 °F); 10 to 90% RH noncondensing	19.7 x 8.9 x 11.4 cm (7.7" x 3.5" x 4.5")	12 VDC	<b>07532-25</b>

\*Intermittent duty: limit output to 18 psig – Continuous duty: limit output to 16 psig.



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.