

# Cole-Parmer®

## Cole-Parmer® See-saw Rockers - SH-200D-S, SH-200D-S-L

- Small space saving design - ideal for personal use
- Gentle See-saw rocking action
- Incubator safe
- Digital speed control and built-in timer
- Optional tier system available to increase capacity



SH-200D-S  
Mini See-saw Rocker



SH-200D-S-L See-saw Rocker

## Cole-Parmer® See-saw Rockers - SH-200D-S, SH-200D-S-L

- Small space saving design or laboratory sized unit
- Digital speed control and built-in timer
- Optional tier system available to increase capacity

The Cole-Parmer® models have a see-saw rocking action that creates a wave motion within vessels such as culture flasks, Petri dishes etc. Digital selection of both speed and time facilitates accurate and reproducible conditions. Both parameters are controlled via an easy to use control dial and displayed on the bright LED display.

Both models are supplied with a non-slip mat which holds flat based vessels securely in place during mixing. An accessory tier system can be fitted in seconds, without the need for tools, tripling the available space for samples. Each tier is securely held in place by magnets.

### Technical Specification

Specification	SH-200D-S	SH-200D-S-L
<b>Rocking action</b>	See-saw wave	See-saw wave
<b>Platform dimensions (w x l)</b>	235 x 235mm	355 x 355mm
<b>Speed range</b>	5 to 70rpm	5 to 70rpm
<b>Angle of tilt</b>	7°	7°
<b>Maximum load</b>	3kg	10kg
<b>Overall dimensions (w x d x h)</b>	240 x 300 x 150	360 x 420 x 170
<b>Tier height, mm</b>	125/tier	125/tier
<b>Operational temp. range</b>	+4 to +40°C	+4 to +40°C
<b>Max. permissible humidity</b>	80%	80%
<b>Net weight</b>	5kg	10kg
<b>Electrical supply</b>	120V, 60 Hz, 50W 230V, 50Hz, 50W	120V, 60 Hz, 50W 230V, 50Hz, 50W

### Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Cole-Parmer® See-saw rocker mini	51901-14	SH-200	SH-200D-S	SSM4
Cole-Parmer® See-saw rocker mini 120V	51901-15	SH-200	SH-200D-S-120	SSM4/120V/60
Cole-Parmer® See-saw rocker large	51901-06	SH-200	SH-200D-S-L	SSL4
Cole-Parmer® See-saw rocker large 120V	51901-07	SH-200	SH-200D-S-L-120	SSL4/120V/60

### Accessories

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Cole-Parmer® Tier system (2 platforms plus 8 bars)	51700-93	SH-200	SH-200D-S/SH200D-S/120	SSM3/1
Cole-Parmer® Tier system (2 platforms plus 8 bars)	51700-94	SH-200	SH-200D-S-L/SH200D-S-L/120	SSL3/1



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.