

# Cole-Parmer®

## Cole-Parmer® Digital Water Baths WB-300-6, WB-300-15, WB-300-24

- Choice of sizes: 6, 15 or 24 Litres
- Digital display and selection of temperature
- Low level water sensor
- Integral drain
- Supplied with robust polycarbonate lid



WB-300-6



WB-300-15



WB-300-24

## Cole-Parmer® Digital Water Baths WB-300-6, WB-300-15, WB-300-24

Three robust and reliable water baths, each with easy to clean stainless steel tank and sturdy metal outer casing. Digital water baths with easy to use temperature control and clear to read LED water temperature. Incorporated over temperature protection system that tracks the set temperature and controls the heater in the event of a fault. The heater is mounted underneath the tank to allow easy cleaning. Includes low level water sensor which cuts power to prevent the bath boiling dry. Incorporated drain for easy emptying of the bath. Bath covers are included as standard with each bath to help reduce heat losses and evaporation. The covers are manufactured from robust polycarbonate allowing full visibility of the bath contents.

### Ordering Information

Description	Ordering Number	Series No.	Model No.	Legacy Sku.
Digital water bath 6L	12122-91	WB-300	WB-300-6	SWB6D
Digital water bath 15L	12122-71	WB-300	WB-300-15	SWB15D
Digital water bath 24L	12122-81	WB-300	WB-300-24	SWB24D
Polycarbonate lid for 6L bath	12122-93	WB-300	N/A	SWB6D/1
Polycarbonate lid for 15L bath	12122-73	WB-300	N/A	SWB15D/1
Polycarbonate lid for 24L bath	12122-83	WB-300	N/A	SWB24D/1
Stainless Steel lid for 6L bath	12122-92	WB-300	N/A	SWB6/1
Stainless Steel lid for 15L bath	12122-72	WB-300	N/A	SWB15/1
Stainless Steel lid for 24L bath	12122-82	WB-300	N/A	SWB24/1

### Technical Specification

Specification	WB-300-6	WB-300-15	WB-300-24
Capacity	6L	15L	24L
Temperature range	25°C to 99.9°C	25°C to 99.9°C	25°C to 99.9°C
Temperature stability	+/- 0.5C	+/- 0.5C	+/- 0.5C
Internal dimensions (w x d x h),mm	300 x 150 x 200	300 x 325 x 200	300 x 500 x 200
Overall dimensions (w x d x h), mm	335 x 230 x 280	335 x 408 x 280	335 x 584 x 280
Net Weight, kg	6.5kg	9.0kg	12.8kg
Heater Power	350W	750W	1400W
Electricity supply	230V, 50Hz,	230V, 50Hz,	230V, 50Hz,
IP Rating	31	31	31

Note: Minimum set temperature on all models is 25°C



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.