


Cole-Parmer® Water Bath Comparison Guide



Specifications	WB-200 Series			WB-300 Series			WB-400 Series				
Item Numbers	78905-20 78905-21	78905-22 78905-23	78905-24 78905-25	12122-91	12122-71	12122-81	12105-84 12105-85	12105-86 12105-87	12105-88 12105-89	12105-90 12105-91	12105-92 12105-93
Capacity	6.6 L	11 L	27 L	6 L	15 L	24 L	2 L	5 L	10 L	20 L	28 L
Temperature Range	Ambient +5 to 100 °C			25 to 99.9 °C			Ambient +5 to 100 °C				
Temperature Stability	±0.2 °C			±0.5 °C			±0.2 °C				
Heater Wattage	600 W	800 W	1000 W	350 W	750 W	1400 W	120 W	360 W	1000 W	1400 W	1400 W
Low Water Sensor	—	—	—	●	●	●	—	—	—	—	—
Drain Valve	●	●	●	●	●	●	●	●	●	●	●
Programmable Timer	●	●	●	—	—	—	●	●	●	●	●
Programmable Presets	—	—	—	—	—	—	●	●	●	●	●
Overtemperature Protection	●	●	●	●	●	●	●	●	●	●	●
Delayed Start	●	●	●	—	—	—	—	—	—	—	—
Visual Alarms	●	●	●	—	—	—	●	●	●	●	●
Multiple Language Interface	—	—	—	—	—	—	●	●	●	●	●
Lid Included	●	●	●	●	●	●	●	●	●	●	●
Electrical	120 VAC, 60 Hz or 220 VAC, 50/60 Hz			230 VAC, 50 Hz			120 VAC, 60 Hz or 240 VAC, 50 Hz				
Warranty	2 years			3 years			2 years				



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.