



Model: CMCF650

650L Laboratory Chest Freezer

Overview

- -10°C to -25°C Temperature Range
- External Digital Control and Temperature Display
- Fitted Handle with Key Lock and 2 Keys
- Twin Baskets to Separate Frozen Items
- 4 Star Freezer Rating

5 YEARS PARTS &
LABOUR WARRANTY

ANNUAL SERVICE &
CALIBRATION AVAILABLE

Key Features



Controller

Controller system allows reading, recording, and resetting with the touch of one button.



Flexible Storage

Multiple baskets to provide a flexible storage solution.



Spark Free

No internal light or fan, suitable for a spark free laboratory.

Details

Model: CMCF650

Gross Capacity (litres)	567
Net Capacity (litres)	560
Unit External Width (mm)	1805
Unit External Depth (mm)	750
Unit External Height (mm)	850
Product Weight (kg)	79
Packed Width (mm)	1850
Packed Depth (mm)	780
Packed Height (mm)	880
Packed Weight (kg)	86
Lid Thickness (mm)	80
Wall Thickness (mm)	70
Temperature Range	-10°C to -25°C
External Finish, Internal Finish	Aluminium, PCM
Battery Backup	No
Door Type, Door Style	Flat, PCM, Hinged
Open Door Sensor	No
Lock	Yes, Integrated with handle
Adjustable Feet (number)	4
Cooling System	Static

Temperature Display on Panel, Temp Alarm, Remote Alarm Contacts, Audible Alarms, Visual Alarms	Yes
Number of Shelves/Baskets	2
Defrost Drain	With
Shelf Material	White Wire
Interior Light	No
Light Type	N/A
Fan	No
Spark Free	Yes
Refrigerant Type	R600a
CFC Free, HCFC Free, Ammonia Free	Yes
Voltage Supply	220-240V/ 50Hz
Certification	CE/UK
Input Power (W)	125
Energy Consumption (KWh/annum)	314
Energy Rating	E

Related Products



CMF86V208
-86°C Ultra Low
Temperature Freezer



CMF86V80
-86°C Ultra Low
Temperature Freezer



CMCF550
550L Laboratory
Chest Freezer



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.