

Model: CMST50

Small Staff Room Fridge

Overview

- +2°C to +10°C temperature range
- Energy Rating A+
- Annual service & Calibration available
- Reversible Door
- Adjustable legs
- 6 Setting temperature control
- Recessed handle
- 2 Wire interior shelf
- Door shelf
- 479W x 447D x 552H mm



**2 YEARS PARTS &
LABOUR WARRANTY**

**ANNUAL SERVICE &
CALIBRATION AVAILABLE**

Key Features



Flexible storage

Multiple shelf locations. Back stop on every shelf provides a flexible storage solution.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.