

# DM0506

## Low-speed Centrifuge



### Available adapters

#### A5P17

Cat. No.19500022

5mL adapter, Adapted to 13×75mm blood collection tube



#### A2P17

Cat. No.19500023

1.5/2mL Cryogenic Vials/ centrifuge tube adapter



## Features

- Speed range from 300-5000 rpm
- Set and displayed by RPM or G-force
- Two programs P1/P2, easy to start the procedure by one key
- Two deceleration speeds setting to meet different requirements
- Easy-to-read LCD display and sound alert
- Automatic lid-lock release after running
- Maintenance-free brushless DC motor

## Applications

DM0506 is widely used in the medical industry. Ideal for separation of serum, plasma, urine, fecal samples in biochemical & medical laboratories. It is an essential equipment of every clinical laboratory, blood bank, medical college and medical research institute.

| Specifications            | DM0506                                                           |
|---------------------------|------------------------------------------------------------------|
| Speed Range               | 300-5000rpm, increment:10rpm                                     |
| Max. RCF                  | 2600×g, increment: 10×g                                          |
| Speed accuracy            | ±20rpm                                                           |
| Rotor type                | 6x15mL/10mL/7mL/1.5-5mL                                          |
| Run Time                  | 30sec-99min, HD (continuous running)                             |
| Display                   | LCD                                                              |
| Preloaded Program         | P1/P2                                                            |
| Acceleration/Deceleration | 1 ↑ 2 ↓                                                          |
| Safety Devices            | Door interlock, Overspeed detector, Automatic internal diagnosis |
| Motor                     | Brushless DC motor                                               |
| Power Requirements        | Single-phase, 100V-240V, 50Hz/60Hz, 70W                          |
| Dimension[WxDxH]          | 300X240X180mm                                                    |
| Weight                    | 5.2kg                                                            |
| Advanced features         | Speed/RCF switch; sound-alert function                           |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.