



Remote Temperature Monitoring Device USP's



- Monitor up to 4 appliances in the same room from one device
- Monitor appliances with a temperature range of minus 40C to plus 100C
- Built-in Ambient temp and Humidity sensors (Measure Temp of the room/environment that the appliances are kept in)
- Device operates on a Global multi-network GSM Sim card(2G-5G)
- Mains-powered with backup battery that lasts for 3 weeks and more in the event of a power outage
- World Health Organisation PQS approved, (Trusted by over 130 International Governments)
- **Connected Communication** - With proactive alerts, you can relax, focus on your work

and be safe in the knowledge that should anything occur with your valuable temperature-controlled contents, you will be notified wherever you may be.

- **Access:** Whether via the desktop portal, or via the app, you can access your appliances' alarm state and related information from wherever you may be.
- **Immediate alarms:** In the event of an alarm status, you will be notified via the app or by email to ensure action can be taken to preserve your contents.
- **Ease of recording:** Temperature data is downloadable in CSV format for your records.

- **Secure:** Temperature data is securely stored on your portal. All data is backed up with failsafe data recovery in place.
- **Simplicity:** You can view at a glance whether your devices are operating as they should from your homepage, or should you need to look a little deeper, there are quick view temperature graphs for every appliance
- **Scalability:** An unlimited number of connected devices can be seen from your portal across any location.
- **Data storage:** There is no need to archive data, you can access all historical data at any time.
- **Accurate history:** You can view the performance history and add comments for any actions taken, ensuring a seamless, recorded data trail.
- **Reliable monitoring:** Appliance temperature is recorded at regular intervals allowing you to view an up-to-date status overview in just a few moments. All reports can be automated via email whether daily, weekly or monthly.

Portal Credentials:

- WHO- PQS Certified (Performance Quality, Safety)
- ISO/IEC27001 - 2022 (Information Security Management System)
- ISO9001 - 2015 (Quality Management System)
- SSL Labs rating: A+



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.