

COUNTER TOP

Glass Door: DLS95PG

Overview

Our intuitive range refrigerators feature built-in wireless monitoring software, enabling you to access the refrigerator's performance remotely.

Plus, our bespoke touchscreen controller, designed specifically for the end user with simple access to key information.



Key Features



Advanced Controller

Easy to use touch screen controller, with no awkward buttons to press



Connected

Wireless temperature monitoring built-in



Dual Temperature Probes

View your daily minimum and maximum temperatures in one touch



Battery Back-Up

A/V alarms should your temperature stray or the door be left open



Antimicrobial Protection

Biocote® handle helping protect from microbes such as E. coli, Campylobacter, Salmonella and MRSA



Self-Closing Doors

Fitted with self-closing hinges to help avoid any temperature loss

DLS95PG - Product Specification

| Counter Top | |
|----------------------------------|--------------------------|
| Model | DLS95PG |
| Gross Capacity (Litres) | 95 |
| Net Capacity (Litres) | 83 |
| Unit External Width (mm) | 520 |
| Unit External Depth (mm) | 630 (excl handle) |
| Unit External Height (mm) | 690 |
| Packed Width (mm) | 530 |
| Packed Depth (mm) | 640 |
| Packed Height (mm) | 690 |
| Packed Weight (kg) | 32 |
| Internal Dimensions (mm) (WxDxH) | 430x390x510 |
| Temperature Range | +2C to +8C |
| External Finish | Pre-coated steel |
| Internal Finish | Plastic & Aluminium |
| Type of Plastic | HIPS (Food grade) |
| Door Type | Right hand Hinged |
| Door Style | Glass / Self Closing |
| Open Door Depth | 1137 |
| Replaceable Door Seal | Yes |
| Open Door Sensor | Yes |
| Lock | Yes / Mechanical |
| Adjustable Feet | One |
| Porthole | Yes |
| Porthole Size | 16mm |
| Temperature Display On Panel | LED / Touchscreen |
| Temperature Alarm | Yes |
| Remote Alarm Contacts | Yes |
| Audible Alarms | Yes |
| Visual Alarms | Yes |
| Number Of Shelves | 2 |
| Number Of Shelf Positions | 3 |
| Shelf Size | 444 x 400 |
| Shelf Material | PVC coated wire |
| Number Of Baskets | 1 |
| Basket Size | 410x250x148 |
| Basket Material | PVC coated wire |
| Interior Light | Yes |
| Light Type | LED |
| Light Activation | Door Switch / Controller |
| Fan | Yes |
| Spark Free | No |
| Defrost Setting | Auto |
| Controller Type | Digital / Wireless |
| Battery Back Up | Yes |
| Controller Connectivity | Wireless / SD Download |
| Number Of Probes | 2 |
| Probe Type | Air / Load |
| Ambient Temperature Range | 10 - 32C |
| Refrigerant Type | R600a |
| Plug Type | UK |
| Voltage Supply | 220-240V |
| Frequency (Hz) | 50Hz |
| Insulation Thickness | 37 |
| Anti Microbial Handle | Yes / BioCote |



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.