

Silenced Compressor Station

187-15 (GBXAG132XTABKM)



- Oil free compressor technology
- Maintenance free except for the intake filter(s) and sinter filter
- Ready for connection
- Adjustable regulator on main control panel
- Internal membrane drier suitable for air purity up to (1.3.1)
- Contained within a temperature controlled acoustic cabinet with lockable door
- Low noise level 58-59 dB (A)
- Suitable for use in an S1 condition up to 30°C ambient
- Medical grade castors for ease of movement
- Energy recycling unit removes need for an external moisture drain
- Extended warranty available (on request)
- Optional regulated distribution manifold for up to 5 applications
- Auxiliary connection to enable external pressure vessel to be added if required



Technical data

Delivery Rate @ 0 bar	105 l/min
Delivery Rate @ 5 bar	62 l/min
Receiver volume	12.5 litres
Voltage	230 V
Frequency	50/60 Hz
Protection class	IP20
Dimensions in cm (H x W x D)	105 x 56 x 45
Weight	80 kg (Unpacked)

:



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.