

REFERENCE : 060976CL

**MP Orbital Shaker for microplates
ClearLine®**

Features :

- Orbital shaker with digital display
- For 4 standard or deepwell microplates
- Suitable for cold rooms and incubator
- Comes with a microplate platform



Vitesse d'agitation (rpm)	200 to 1500
Incrément (rpm)	10
Orbite (mm)	3
Capacity	4 microplates
Shaker Dimensions L x P x H (mm)	330 x 300 x 120
Affichage	Digital
Operating Temperature	4 à 65°C
weight (kg)	5,5
Electrical power	230V, 50-60Hz

Sales Unit :

1 shaker



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.