

Logiclave

Autoclave & Steam Sterilizer Technology



LogiclaveLAB Standard Product Range

Model Number	Volume (L)	Double Door Option	Chamber shape	Door Movement	Chamber Dimention H x W x D (mm)	Power Supply
Logiclave LAB100V	100	No	Cylindrical	Vertical	500 x 500 x 500	3 Phase, 400VAC, 3-4 kW
Logiclave LAB145V	145	No	Cylindrical	Vertical	500 x 500 x 750	3 Phase, 400VAC, 3-4 kW
Logiclave LAB150V	150	No	Rectangular	Vertical	600 x 500 x 500	3 Phase, 400VAC, 3-4 kW
Logiclave LAB220V	220	No	Cylindrical	Vertical	500 x 500 x 1000	3 Phase, 400VAC, 3-4 kW
Logiclave LAB225V	225	No	Rectangular	Vertical	600 x 500 x 750	3 Phase, 400VAC, 3-4 kW
Logiclave LAB300V	300	Yes	Rectangular	Vertical	600 x 500 x 1000	3 Phase, 400VAC, 3-4 kW
Logiclave LAB375V	375	Yes	Rectangular	Vertical	600 x 500 x 1250	3 Phase, 400VAC, 3-4 kW
Logiclave LAB435V	435	Yes	Rectangular	Vertical	660 x 660 x 1000	3 Phase, 400VAC, 3-4 kW
Logiclave LAB585V	585	Yes	Rectangular	Vertical	660 x 660 x 1350	3 Phase, 400VAC, 3-4 kW
Logiclave LAB600	600	Yes	Rectangular	Horizontal	920 x 660 x 1000	3 Phase, 400VAC, 3-4 kW
Logiclave LAB650V	650	Yes	Rectangular	Vertical	660 x 660 x 1500	3 Phase, 400VAC, 3-4 kW
Logiclave LAB800	800	Yes	Rectangular	Horizontal	920 x 660 x 1350	3 Phase, 400VAC, 3-4 kW
Logiclave LAB900	900	Yes	Rectangular	Horizontal	920 x 660 x 1500	3 Phase, 400VAC, 3-4 kW
Logiclave LAB960	960	Yes	Rectangular	Horizontal	1200 x 800 x 1000	3 Phase, 400VAC, 3-4 kW
Logiclave LAB1295	1295	Yes	Rectangular	Horizontal	1200 x 800 x 1350	3 Phase, 400VAC, 3-4 kW
Logiclave LAB1440	1440	Yes	Rectangular	Horizontal	1200 x 800 x 1500	3 Phase, 400VAC, 3-4kW
Logiclave LAB1900	1900	Yes	Rectangular	Horizontal	1200 x 800 x 2000	3 Phase, 400VAC, 3-4kW
Logiclave LAB2565	2565	Yes	Rectangular	Horizontal	1500 x 950 x 1800	3 Phase, 400VAC, 4-6kW
Logiclave LAB2850	2850	Yes	Rectangular	Horizontal	1500 x 950 x 2000	3 Phase, 400VAC, 4-6kW
Logiclave LAB3135	3135	Yes	Rectangular	Horizontal	1500 x 950 x 2200	3 Phase, 400VAC, 4-6kW
Logiclave LAB3420	3420	Yes	Rectangular	Horizontal	1500 x 950 x 2400	3 Phase, 400VAC, 4-6kW
Logiclave LAB3590	3590	Yes	Rectangular	Horizontal	2100 x 950 x 1500	3 Phase, 400VAC, 4-6kW
Logiclave LAB3990	3990	Yes	Rectangular	Horizontal	2100 x 950 2000	3 Phase, 400VAC, 4-6kW
Logiclave LAB4390	4390	Yes	Rectangular	Horizontal	2100 x 950 x 2200	3 Phase, 400VAC, 4-6kW
Logiclave LAB4725	4725	Yes	Rectangular	Horizontal	2100 x 1250 x 1800	3 Phase, 400VAC, 4-6kW
Logiclave LAB4780	4780	Yes	Rectangular	Horizontal	2100 x 950 x 2400	3 Phase, 400VAC, 4-6kW
Logiclave LAB5250	5250	Yes	Rectangular	Horizontal	2100 x 1250 x 2000	3 Phase, 400VAC, 4-6kW
Logiclave LAB5775	5775	Yes	Rectangular	Horizontal	2100 x 1250 x 2200	3 Phase, 400VAC, 4-6kW
Logiclave LAB6300	6300	Yes	Rectangular	Horizontal	2100 x 1250 x 2400	3 Phase, 400VAC, 4-6kW



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.