



## VARIABLE VOLUME MECHANICAL PIPETTES



### I Ultimate Performance

Built for traceability and comfort. PIPETMAN L is designed to increase daily productivity, giving you the ability to pipette longer without fatigue.

### I Secure Pipetting

Reduce pipetting errors with volume-locking. Avoid accidental volume changes during pipetting cycles, ensuring your accuracy and precision.

### I Optimized Comfort

Lightweight and light pipetting forces. Right-handed or left-handed, no problem. Adjust the location of the tip ejector to the most comfortable position for you.





## Classic metal or plastic tip ejector

Choose between a stainless steel tip ejector for more robustness or a plastic tip ejector to avoid corrosion.



## Choose your model

Choose your channel number and volume range. Each ensures consistent volume levels, allowing you to aspirate and dispense your samples with confidence.

PIPETMAN L Single Channel — Variable Volume Models				
Model	PIPETMAN DIAMOND Tips		Part Number	Volume (μL)
<b>P2L</b>	D10 DL10	DF10 DFL10	<b>FA10001P</b> <b>FA10001M</b>	0.2 - 2
<b>P10L*</b>	D10 DL10	DF10 DFL10	<b>FA10002P</b> <b>FA10002M</b>	1 - 10
<b>P20L</b>	D200	DF30	<b>FA10003P</b> <b>FA10003M</b>	2 - 20
<b>P100L</b>	D200	DF100	<b>FA10004P</b> <b>FA10004M</b>	10 - 100
<b>P200L</b>	D200 D300	DF200 DF300	<b>FA10005P</b> <b>FA10005M</b>	20 - 200
<b>P1000L</b>	D1000 D1200	DF1000 DF1200	<b>FA10006P</b> <b>FA10006M</b>	100 - 1000
<b>P5000L</b>	D5000		<b>FA10007**</b>	500 - 5000
<b>P10mLL</b>	D10mL		<b>FA10008**</b>	1000 - 10000
PIPETMAN L Multichannel				
<b>P8x10L</b>	D10 DL10 DF10 DFL10		<b>FA10013</b>	0.5 - 10
<b>P12x10L</b>			<b>FA10014</b>	
<b>P8x20L</b>	DL10 DFL10 DF30 D200		<b>FA10009</b>	2 - 20
<b>P12x20L</b>			<b>FA10010</b>	
<b>P8x200L</b>	D200 D300 DF200 DF300		<b>FA10011</b>	20 - 200
<b>P12x200L</b>			<b>FA10012</b>	
<b>P8x300L</b>	D200 D300 DF200 DF300		<b>FA10015</b>	20 - 300
<b>P12x300L</b>			<b>FA10016</b>	

For a pipette with a plastic tip ejector, the part number ends with the letter P. For a pipette with a stainless steel tip ejector, the part number ends with the letter M.



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.