



## Technical information & Ordering data

### 96-well

#### Standard microplates



| Bottom                          | U-bottom                                 |                                        | V-bottom                                 |                                        | F-bottom                                 |                                        | C-bottom                                 |                                        |
|---------------------------------|------------------------------------------|----------------------------------------|------------------------------------------|----------------------------------------|------------------------------------------|----------------------------------------|------------------------------------------|----------------------------------------|
| Color                           | transparent                              |                                        | transparent                              |                                        | transparent                              |                                        | transparent                              |                                        |
| Well volume [μl]                | 330                                      |                                        | 360                                      |                                        | 350                                      |                                        | 350                                      |                                        |
| Working volume [μl]             | 40-300                                   |                                        | 40-330                                   |                                        | 50-320                                   |                                        | 50-33                                    |                                        |
| Bottom thickness [μm]           | 850                                      |                                        | 850                                      |                                        | 850                                      |                                        | 850                                      |                                        |
| Well surface [mm <sup>2</sup> ] | n/a                                      |                                        | 33                                       |                                        | 32                                       |                                        | 25                                       |                                        |
|                                 |                                          | sterile                                |                                          | sterile                                |                                          | sterile                                |                                          | sterile                                |
| Lid                             | 20 pieces<br>(1 lid/stack)               | 50 pieces                              | 20 pieces<br>(1 lid/stack)               | 50 pieces                              | 20 pieces<br>(1 lid/stack)               | 50 pieces                              | 20 pieces<br>(1 lid/stack)               | 50 pieces                              |
| Pack of                         | 100 pieces<br>(20 stacks<br>of 5 plates) | 50 pieces<br>(individually<br>wrapped) | 100 pieces<br>(20 stacks<br>of 5 plates) | 50 pieces<br>(individually<br>wrapped) | 100 pieces<br>(20 stacks<br>of 5 plates) | 50 pieces<br>(individually<br>wrapped) | 100 pieces<br>(20 stacks<br>of 5 plates) | 50 pieces<br>(individually<br>wrapped) |
| Cat. No.                        | 781600                                   | 781660                                 | 781601                                   | 781661                                 | 781602                                   | 781662                                 | 781603                                   | 781663                                 |



| Bottom                          | U-bottom                              |  | F-bottom                              |                                     | U-bottom                              |  | F-bottom                              |                                     |
|---------------------------------|---------------------------------------|--|---------------------------------------|-------------------------------------|---------------------------------------|--|---------------------------------------|-------------------------------------|
| Color                           | white                                 |  | white                                 |                                     | black                                 |  | black                                 |                                     |
| Well volume [μl]                | 330 μ                                 |  | 350                                   |                                     | 330                                   |  | 350                                   |                                     |
| Working volume [μl]             | 40-300                                |  | 50-320                                |                                     | 40-300                                |  | 50-320                                |                                     |
| Bottom thickness [μm]           | 850                                   |  | 850                                   |                                     | 850                                   |  | 850                                   |                                     |
| Well surface [mm <sup>2</sup> ] | n/a                                   |  | 32                                    |                                     | n/a                                   |  | 32                                    |                                     |
|                                 |                                       |  |                                       | sterile                             |                                       |  |                                       | sterile                             |
| Lid                             | 20 pieces<br>(1 lid/stack)            |  | 20 pieces<br>(1 lid/stack)            | 50 pieces                           | 20 pieces<br>(1 lid/stack)            |  | 20 pieces<br>(1 lid/stack)            | 50 pieces                           |
| Pack of                         | 100 pieces<br>(20 stacks of 5 plates) |  | 100 pieces<br>(20 stacks of 5 plates) | 50 pieces<br>(individually wrapped) | 100 pieces<br>(20 stacks of 5 plates) |  | 100 pieces<br>(20 stacks of 5 plates) | 50 pieces<br>(individually wrapped) |
| Cat. No.                        | 781604                                |  | 781605                                | 781665                              | 781607                                |  | 781608                                | 781668                              |

Lid for 96-well standard plate, see page 84:  
without condensation rings Cat. No. 782151  
with condensation rings Cat. No. 782150



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.