


MMS-3000 mini magnetic stirrer

Compact magnetic stirrer for routine laboratory procedures.

- Adjustable speed control: 0 to 3000rpm
- Maximum stirring volume: 20 litres



- Applications:**
- PH metering, extraction and dialysing with small quantities of substances
 - Any requirement for stirring solutions up to 20 litres

Stirrers – models and specifications

MMS-3000

Mini magnetic stirrer



h: 75mm
d: 230mm
w: 185mm
weight: 1.5kg

Operating temperature range	°C	+4 to 40
Speed	rpm	0 to 3000
Liquid stirring viscosity	mPa	Up to 1170
Stirring volume	L	20
Working plate size	mm	ø 160
External power supply		Input AC 120-230V, 50/60Hz Output DC 12V
Input voltage	V dc	12
Power consumption	W	3 (0.3A)



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.