


CENTRIFUGE SELECTION GUIDE

MODEL	MAX. RPM	RCF	CONTAINERS																										
			8-tube MICRO-TUBE PER STRIP	0.2-2.0 ml MICRO-TUBE	1.5 ml CRYO-TUBE	2.0 ml MICROFILTER-TUBE	5.0 ml MICRO-TUBE	3-5 ml TUBE & VACUTAINER	7 ml TUBE & VACUTAINER	10 ml TUBE & VACUTAINER	15 ml TUBE	15 ml conical TUBE	50 ml TUBE	50 ml conical TUBE	50 ml conical TUBE	95 ml TUBE	100 ml TUBE	250 ml TUBE	250 ml (C) BOTTLE	500 ml BOTTLE	500 ml (C) BOTTLE	750 ml BOTTLE	1,000 ml BOTTLE	MICROPLATE	DEEP-WEEL PLATE	FILTER PLATE	10-20 ml SYRINGE	60 ml SYRINGE	
mini 6	6,500	2,431	2	6	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
mini	13,500	12,300	4	12	10	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
1536	15,000	21,583	4	36	-	-	10	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
MICRO	17,000	27,237	8	36	-	18	12	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
1730R	13,000	16,343	-	16	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
1848R	18,000	30,536	8	48	-	-	20	-	-	-	-	-	-	12	4	8	4	-	-	-	-	-	-	-	-	-	-	-	
406	4,000	2,075	-	-	-	-	-	12	12	12	12	6	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
416	4,000	2,700	-	-	-	-	-	16	16	16	16	16	10	10	10	-	-	-	-	-	-	-	-	-	-	-	-	-	
624R	4,800	3,888	-	24	-	-	-	32	32	32	32	12	8	8	4	4	4	-	-	-	-	-	-	-	-	-	-	-	
1248/1248R	12,000	16,582	8	30	-	-	-	24	24	24	24	24	6	8	-	6	-	-	-	-	-	-	-	-	-	-	-	-	
1580	5,000	5,394	-	80	-	-	-	48	48	48	48	36	16	12	16	8	6	4	4	-	-	-	-	-	16	4	-	8	4
1580R	15,000	25,910	8	30	-	-	-	24	24	24	24	24	8	8	-	6	-	6	-	-	-	-	-	-	-	-	-	-	-
1696R	16,000	27,590	8	48	-	-	20	-	-	8	24	12	8	8	-	6	-	6	-	-	-	-	-	-	-	-	-	-	-
FLOOR-STANDING, HIGH-SPEED	12,000	16,582	8	30	-	-	-	24	24	24	24	24	8	8	-	6	-	-	-	-	-	-	-	-	-	-	-	-	-
1236R	5,000	5,394	-	80	-	-	-	48	48	48	48	36	16	12	16	8	6	4	4	-	-	-	-	-	16	4	-	8	4
1736R	17,000	32,310	8	30	-	-	-	12	12	12	8	8	-	6	-	6	6	6	6	-	6	-	-	-	-	-	-	-	-
2236R	22,000	54,163	8	30	-	-	-	12	12	12	8	8	-	6	-	6	6	6	6	-	6	-	-	-	-	-	-	-	-



EXPERTS IN LABORATORY EQUIPMENT

- CENTRIFUGE
- CONCENTRATOR
- FREEZE DRYER
- DEEP FREEZER
- CO₂ INCUBATOR
- MIXERS & SHAKER
- CONSUMABLES

GYROZEN PRODUCT OVERVIEW
CENTRIFUGE





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.