

Combined Refrigerator and Freezer

HYCD-282



Refrigerator & Freezer Combination

Independent display for refrigerator section and freezer section



Alarm Functions

High/low temperature, power failure, sensor error, low battery, high/low ambient temperature, and remote alarm



Standard Stainless Steel Interior

Both refrigerator and freezer chamber are stainless steel interior, durable, easy to clean and anti-corrosion

Ergonomic Design



Lockable design



Intelligent alarm



Foot/Caster design



Test hole design

Specifications

HYCD-282

Technical Data

Type	Combined Refrigerator and Freezer	
Cabinet Type	Upright	
Climate Class	4	
Cooling Type	Refrigerator: Forced air cooling Freezer: Direct cooling	
Defrost Mode	Refrigerator: Auto Freezer: Manual	
Refrigerant	CFC-Free	
Noise Level (dB(A))	43	

Performance

Temperature Range(°C)	Refrigerator: 2~8 Freezer: -20~-40
-----------------------	---------------------------------------

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220~240/50	
Power(W)	400	
Electrical Current(A)	2.8	

Dimensions

Capacity(L/Cu.Ft)	Refrigerator:185/6.5 Freezer:97/3.4	
Net/Gross Weight(approx)	145/160(kg)	319.7/352.7(lbs)
Interior Dimensions(W*D*H)	Refrigerator:605*510*720(mm) Freezer:515*465*440(mm)	23.8*20.0*28.3(in) 20.3*18.3*17.3(in)
Exterior Dimensions(W*D*H)	736*660*1810(mm)	29.0*26.0*71.3(in)
Packing Dimensions(W*D*H)	775*720*1930(mm)	30.5*28.3*76(in)
Container load (20'/40'/40'H)	21/45/45	



Functions

High/Low Temperature	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves/Drawers	Refrigerator: 3/1 Freezer: 0/2
USB Interface	Y

Others

Certificate	CE
-------------	----



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.