

Haier Biomedical

Intelligent Protection of Life Science

Pharmacy Refrigerators

HYC-118A



Air Cooling, Auto Defrost

Better cabinet temperature
uniformity, auto defrost



Five Alarm Functions

High/low temperature,
power failure, door ajar,
sensor error, low battery



Audible and Visual Alarms

Audible buzzer and visible
flashing light

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Specifications

HYC-118A

Technical Data

Cabinet Type	Upright,glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	CFC-Free
Noise((dB(A))	32

Performance

Temp Range(℃)	2~8
---------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50	115/60
Power(W)	165	140
Electrical Current(A)	1.16	1.46

Dimensions

Capacity(L/Cu.Ft)	118/4.2
Net/Gross Weight(approx)	46/51(kg) 101.4/112.4(lbs)
Interior Dimensions(W*D*H)	515*415*630(mm) 20.3*16.3*24.8(in)
Exterior Dimensions(W*D*H)	597*635*835(mm) 23.5*25.0*32.9(in)
Packing Dimensions(W*D*H)	680*690*910(mm) 26.8*27.2*35.8(in)
Container Load(20'/40'/40'H)	48/102/102



Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote System	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves /Basket	3/1

Others

Certificate	CE,UL
-------------	-------



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.