



## Eco smartFunnel™

### Paper Funnel

A clean transfer every time

- Supplied flat-pack to minimize storage
- With a minimal amount of finger pressure the funnel takes its working shape
- Maximum flexibility – cut the stem to any size or shape to meet with your needs
- Fits a wide variety of openings and can be used with virtually any liquid\* or powder

\* Not recommended for hot liquids

Pack of 25

Item No.	Capacity mL	Top Dia. in <sup>2</sup>	Top Dia. cm <sup>2</sup>	Bottom Dia. in <sup>2</sup>	Bottom Dia. cm <sup>2</sup>	Height in	Height cm
120219	186	2.4	6	0.3	0.6	5	12.8

Eco smartFunnel™ is a trademark of LevGo Inc.



## Lead Rings

Stabilize vessels and labware on the bench or in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit

Item No.	Shape	Color	Flask Sizes mL	ID in	ID cm	OD in	OD cm	Weight g
HS8881A	C	Blue ■	125–500	1.9	4.2	3	7.5	215
HS8882A	O			1.9	4.8	3	7.7	210
HS8881B	C	Blue ■	250–1000	2	5.1	3.5	9	465
HS8882B	O			2	5.1	3.5	9	490
HS8881C	C	Blue ■	500–2000	2.1	5.4	3.9	9.8	595
HS8882C	O			2.2	5.7	3.9	9.8	610
HS8881D	C	Blue ■	1000–4000	2.6	6.6	4.4	11.3	890
HS8882D	O			2.8	7	4.5	11.5	935

\*Weights may vary by +/- 3%



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.