



Sprout® Plus Mini Centrifuge

Mini but Mighty - ideal for small runs

Easy to Use

- Tool-free SnapSpin rotors and quick-spin options make the centrifuges fast and easy to use

Ergonomic Design

- The lid was functionally designed to fit in the palm of the hand for one handed lock and unlock allowing for quick tube extraction
- Large lid size accommodates longer and skirted tubes and tall tube lids
- The housing is made from ABS plastic and lid is polycarbonate plastic enabling the unit to be easily cleaned

Safety Measures

- Fitted with a braking function, the unit quickly stops spinning when powered off or when lid is opened
- A weighted, rubberized base eliminates potential movement and dampens vibration and noise



Accessory Kit Included



Specifications

Speed Range	Approx. 6,000 rpm (fully loaded)
RCF	2,000 x g



Sprout® Plus

Item No.	Color	L x W x H in	L x W x H cm
120610	Green	6 x 5 x 4.3	15.3 x 12.8 x 10.9
120611	Blue		

Replacement parts

HS100011	Sprout® Plus Rotor SnapSpin 6 x 1.5/2.0 mL Tube
HS100012	Sprout® Plus Rotor SnapSpin 2 x 8 Strip, 0.2 mL
120148	Sprout® Plus Lid O-Ring, Pack of 10
120491	Sprout® Plus Power Adapter
120015	Tube Adapters, 0.2 mL, 0.5 mL, Pack of 12 (6 ea.)

Package Includes

Sprout® Plus Mini Centrifuge
6-Place 1.5/2.0 mL Tube Rotor
16-Place 0.2 mL Tube Rotor (singles, or 2 x 8 strip)
6 Tube Adapters 0.2 mL
6 Tube Adapters 0.5 mL
4 Lid O-Rings
Storage Case for Rotors and Adapters
Tube Rack
1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



Need tubes?

See our portfolio of 0.2 mL, 0.5 mL, 1.5 mL Microcentrifuge Tubes





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.