



## Cardboard Slide Tray

Cardboard slide trays hold standard slides (76 x 26 mm and 3 x 1 in)

Cardboard slide folders offer a convenient solution for safely working, transporting and short-term storing of slides

### Organize & Protect

- Index on the front cover corresponds with the number grid inside for identification
- Recessed slide compartments and inside covers for specimen protection
- Slides cannot be accessed when cover is closed reducing contamination

### Sturdy Construction

- Compartment holder is made from a double layer of densely pressed cardboard
- Cover is single layer cardboard
- Front protection flaps are attached with durable tape



Item No.	Color	Capacity	L x W x H in	L x W x H cm
HS9902	Black			
HS9913	Blue			
HS9914	Green	20-Place	8.2 x 13.4 x 0.4	20.8 x 34.1 x 0.9
120196	Red			
120197	Yellow			
HS9930	Black	30-Place	8.2 x 12.9 x 0.4	28.1 x 32.8 x 1

## Slide Tray Clasps

Effortlessly locate, retrieve and keep slide folders organized

- Assorted color clasps are ideal for color coding identification and to keep slide folder flaps securely closed
- Clasps can be labeled with a permanent marker or a label can be attached

### Slide Tray Clasps

Pack of 100

Item No.	Color	L x W x H in	L x W x H cm
120653	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	2.2 x 1.7 x 0.6	5.6 x 4.4 x 1.6
120654	Clear <input type="checkbox"/>		

*\*Assorted colors (Clear, Grey, Blue, Green, Red)*





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.