



Ergonomic

Slight downward slant to tip ejector button reduces hand fatigue, incidence of repetitive strain injury (RSI)

Volume Adjustment
and low-force dispensing
all in one button for
convenience

**Finger
hook**

**3-Position
volume display**

Volume Locking Ring
prevents accidental
volume changes

Color Coded & Labeled
with volume range (Both top & side
colored components for quick visual
confirmation)

Precision Metal Components
for both shaft and barrel ensure reliable
performance and long life

Pearl® Adjustable Volume Pipettes


Precision liquid handling

A complete range of single- and multi-channel pipette volumes to cover all your possible needs

The sleek, ergonomic design and reduced plunger pressure help reduce the chance of Repetitive Stress Injury (RSI), as the pipettes require 2/3 the operation pressure of most standard pipettes

Each pipette is calibrated and a calibration certificate is supplied. Manufactured in accordance with ISO8655 and GMP regulations

- Internal metal parts for increased accuracy and calibration stability
- Separate tip ejector with the exception of the 500 mL (120105) and 1,000 mL (120106) volume units
- Large digit volume display
- Manufactured in accordance with ISO8655 and GMP regulations



CE

No. of Channels	Item No.	Color	Volume Range μL
Single	120097	Green	0.1–2 μL
	120098	Lt. Blue	0.5–10 μL
	120099	Lt. Purple	2–20 μL
	120100	Orange	5–50 μL
	120101	Blue	10–100 μL
	120102	Yellow	20–200 μL
	120103	Purple	30–300 μL
	120104	Blue	100–1,000 μL
	120105	Pink	500–5,000 μL
	120106	Red	1,000–10,000 μL
Eight	120215	Orange	5–50 μL
	120216	Purple	30–300 μL
Twelve	120217	Orange	5–50 μL
	120218	Purple	30–300 μL



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.