



XT Series

**Extended Time Cryogenic
Small Freezers**

3XTL

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	3XTL
LN ₂ Capacity (L)	3
Static Hold Time (Days)	27
Evap. Rate (Liters/Day)	0.11
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	7.6 / 193
Height (in/mm)	17.2 / 436.8
WEIGHT	
Empty (lb/kg)	7.2 / 3.2
Full Weight (lb/kg)	12.5 / 5.7
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A

8XTL

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	8XTL
LN ₂ Capacity (L)	8
Static Hold Time (Days)	80
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	15.6 / 396.2
Height (in/mm)	19 / 482.6
WEIGHT	
Empty (lb/kg)	19.6 / 8.9
Full Weight (lb/kg)	33.8 / 15.3
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A

10XT

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



SPECIFICATIONS	10XT
LN ₂ Capacity (L)	10
Static Hold Time (Days)	100
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	11.4 / 289.6
Height (in/mm)	23.8 / 604.5
WEIGHT	
Empty (lb/kg)	16.5 / 7.5
Full Weight (lb/kg)	34.3 / 15.6
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150

FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid



20XT

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	20XT
LN ₂ Capacity (L)	20.7
Static Hold Time (Days)	230
Evap. Rate (Liters/Day)	0.09
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	15.6 / 396.2
Height (in/mm)	25.8 / 655.3
WEIGHT	
Empty (lb/kg)	26 / 11.8
Full Weight (lb/kg)	62.8 / 28.5
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150

34XT

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	34XT
LN ₂ Capacity (L)	34
Static Hold Time (Days)	340
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	18.8 / 477.5
Height (in/mm)	26.3 / 668
WEIGHT	
Empty (lb/kg)	34.75 / 15.8
Full Weight (lb/kg)	95.3 / 43.2
CAPACITY	
Number of Canisters	6
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	150



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.