



designed for scientists



## C-MAG MS 4

/// Data Sheet

Magnetic stirrer without heating. With ceramic set-up plate which offers excellent chemical resistance.

- Powerful motor for stirring quantities of up to 5 l (H<sub>2</sub>O)
- Elevated control panel for protection against leaking liquids
- Speed range: 100 - 1500 rpm



designed for scientists

### Technical Data

|                                                                     |                 |
|---------------------------------------------------------------------|-----------------|
| Number of stirring positions                                        | 1               |
| Stirring quantity max. per stirring position (H <sub>2</sub> O) [l] | 5               |
| Motor rating output [W]                                             | 1.5             |
| Direction of rotation                                               | left            |
| Speed control                                                       | scale 0 - 6     |
| Speed range [rpm]                                                   | 100 - 1500      |
| Set-up plate material                                               | ceramic         |
| Set-up plate dimensions [mm]                                        | 100 x 100       |
| Dimensions (W x H x D) [mm]                                         | 150 x 105 x 260 |
| Weight [kg]                                                         | 2.716           |
| Permissible ambient temperature [°C]                                | 5 - 40          |
| Permissible relative humidity [%]                                   | 80              |
| Protection class according to DIN EN 60529                          | IP 21           |
| Voltage [V]                                                         | 230 / 120 / 100 |
| Frequency [Hz]                                                      | 50/60           |
| Power input [W]                                                     | 30              |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.