


96 Well PCR Plate



PRODUCT SPECIFICATIONS

Catalog Numbers:	IST-401-096MP IST-401-096TP	Pk of 100 Pk of 10	NON-STERILE NON-STERILE
Product Description:	0.2mL 96-well PCR plate Non-skirted, Alpha-numeric labeling, natural-color		
Product Dimensions:	Length: 117.856mm (4.64") Width: 77.724mm (3.06") Height: 20.828mm (.82")		
Packaging:	10 plates per pack 10 packs per case (100 plates total)		
MicroBiological testing:	Certified testing is performed lot by lot. These tests include Nuclease(Rnase Dnase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.		
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.		



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.