

Jeio Tech 3.5 Litre Stirred Water Bath (BW-05H)

Product Code: AAH45115K

Brand: Jeio Tech

The Jeio Tech 3.5 Litre Stirred Water Bath (BW-05H) features a microprocessor PID control with automatic tuning and calibration plus digital LED displays with a quick adjustable dial knob. The powerful circulation system, using the patented agitator, ensures temperature accuracy of $\pm 0.1^{\circ}\text{C}$.

Corrosion resistant stainless steel construction allows the use of other bath fluids as well as water and silicone oil. Side handles are included for easy transport. Safety is ensured with over-temperature and over-current protection. A keypad lock prevents accidental changes.

Features

Performance

- Temperature range from ambient $+5^{\circ}\text{C}$ to 100°C
- Microprocessor PID control/automatic tuning/calibration
- Powerful circulation system using our patented agitator

Convenience

- Corrosion resistant stainless steel construction allowing other bath fluids as well as water and silicone fluids
- Digital LED displays with quick adjustable dial knob (0.1°C resolution)
- Low-profile design with sloped control panel for easier access and better safety. Also side handles for easy transport

Safety

- Isolated design of heater and sensor by a baffle plate to protect users and samples against direct contact
- Complete safety protection system with warning alarms. Over-temperature, over-current protection
- Keypad lock function prevents accidental changes



Specifications

Bath Volume (L)	3.5
Temperature Range (°C)	Amb. +5 to 100
Temperature Accuracy (°C)	± 0.1
Bath Opening/Depth (W x L, D)	240 x 136, 150mm
Dimensions (W x D x H)	307 x 216 x 266mm
Weight (kg)	6.5



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.