

## Jeio Tech 20 Litre Stirred Water Bath (BW-20H)

**Product Code:** AAH45135K

**Brand:** Jeio Tech

The Jeio Tech 20 Litre Stirred Water Bath (BW-20H) features a microprocessor PID control with automatic tuning and calibration plus digital LED displays with a quick adjustable dial knob. The powerful circulation system, using the patented agitator, ensures temperature accuracy of  $\pm 0.1^{\circ}\text{C}$ .

Corrosion resistant stainless steel construction allows the use of other bath fluids as well as water and silicone oil. Side handles are included for easy transport. Safety is ensured with over-temperature and over-current protection. A keypad lock prevents accidental changes.

### Features

#### Performance

- Temperature range from ambient  $+5^{\circ}\text{C}$  to  $100^{\circ}\text{C}$
- Microprocessor PID control/automatic tuning/calibration
- Powerful circulation system using our patented agitator

#### Convenience

- Corrosion resistant stainless steel construction allowing other bath fluids as well as water and silicone fluids
- Digital LED displays with quick adjustable dial knob ( $0.1^{\circ}\text{C}$  resolution)
- Low-profile design with sloped control panel for easier access and better safety. Also side handles for easy transport

#### Safety

- Isolated design of heater and sensor by a baffle plate to protect users and samples against direct contact
- Complete safety protection system with warning alarms. Over-temperature, over-current protection
- Keypad lock function prevents accidental changes



# Specifications

|                               |                   |
|-------------------------------|-------------------|
| Bath Volume (L)               | 20                |
| Temperature Range (°C)        | Amb. +5 to 100    |
| Temperature Accuracy (°C)     | ± 0.1             |
| Bath Opening/Depth (W x L, D) | 498 x 300, 200mm  |
| Dimensions (W x D x H)        | 564 x 372 x 318mm |
| Weight (kg)                   | 15.5              |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.