

TECHNICAL FEATURE

7NK 91NHDP 01 905



The  :

- Ergonomy and design
- Eliminates the fatigue
- Polyvalent
- Auto-adaptive
- Ecologic

KONCEPT SEATS / POLYURETHANE

Black steel, 5-branch, scuff-resistant, reinforced base with articulated glides

Black gas-lift column with height adjustment from 610 to 855 mm

and shock absorber

Height adjustable aluminium foot-ring

Deformable seat in black KANGO soft polyurethane, tilted of 2° forward

Assembly on silentbloc. Auto-adaptive system in the movements of the user

Lumbar support in black KANGO polyurethane, adjustable in depth



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.