

## Test weights class F2

### Class F2 Milligram weights, flat polygonal sheet, aluminium/German silver

Test weight material: Aluminium 1 mg – 5 mg/German silver 10 mg – 500 mg

Container material: Lined plastic



| Milligram weight |        |          | + | Container   |  | + | DAkkS certificate |  | = | Package price |
|------------------|--------|----------|---|-------------|--|---|-------------------|--|---|---------------|
| KERN             |        | Tol ± mg |   | KERN        |  |   | KERN              |  |   | KERN          |
| 338-01           | 1 mg   | 0,06     |   | 347-009-400 |  |   | 962-451           |  |   |               |
| 338-02           | 2 mg   | 0,06     |   | 347-009-400 |  |   | 962-452           |  |   |               |
| 338-03           | 5 mg   | 0,06     |   | 347-009-400 |  |   | 962-453           |  |   |               |
| 338-04           | 10 mg  | 0,08     |   | 347-009-400 |  |   | 962-454           |  |   |               |
| 338-05           | 20 mg  | 0,10     |   | 347-009-400 |  |   | 962-455           |  |   |               |
| 338-06           | 50 mg  | 0,12     |   | 347-009-400 |  |   | 962-456           |  |   |               |
| 338-07           | 100 mg | 0,16     |   | 347-009-400 |  |   | 962-457           |  |   |               |
| 338-08           | 200 mg | 0,20     |   | 347-009-400 |  |   | 962-458           |  |   |               |
| 338-09           | 500 mg | 0,25     |   | 347-009-400 |  |   | 962-459           |  |   |               |

### Class F2 Individual weights, cylindrical shape, finely turned stainless steel

Test weight material: finely turned stainless steel

Container material: Lined plastic or wooden box (337-150-200 and 337-160-200)



| Individual weight |       |          | + | Container   |  | + | DAkkS certificate |  | = | Package price |
|-------------------|-------|----------|---|-------------|--|---|-------------------|--|---|---------------|
| KERN              |       | Tol ± mg |   | KERN        |  |   | KERN              |  |   | KERN          |
| 337-01            | 1 g   | 0,3      |   | 347-030-400 |  |   | 962-431           |  |   |               |
| 337-02            | 2 g   | 0,4      |   | 347-030-400 |  |   | 962-432           |  |   |               |
| 337-03            | 5 g   | 0,5      |   | 347-030-400 |  |   | 962-433           |  |   |               |
| 337-04            | 10 g  | 0,6      |   | 347-050-400 |  |   | 962-434           |  |   |               |
| 337-05            | 20 g  | 0,8      |   | 347-050-400 |  |   | 962-435           |  |   |               |
| 337-06            | 50 g  | 1,0      |   | 347-070-400 |  |   | 962-436           |  |   |               |
| 337-07            | 100 g | 1,6      |   | 347-070-400 |  |   | 962-437           |  |   |               |
| 337-08            | 200 g | 3,0      |   | 347-080-400 |  |   | 962-438           |  |   |               |
| 337-09            | 500 g | 8,0      |   | 347-090-400 |  |   | 962-439           |  |   |               |
| 337-11            | 1 kg  | 16       |   | 347-110-400 |  |   | 962-441           |  |   |               |
| 337-12            | 2 kg  | 30       |   | 347-120-400 |  |   | 962-442           |  |   |               |
| 337-13            | 5 kg  | 80       |   | 347-130-400 |  |   | 962-443           |  |   |               |
| 337-14            | 10 kg | 160      |   | 347-140-400 |  |   | 962-444           |  |   |               |
| 337-15            | 20 kg | 300      |   | 337-150-200 |  |   | 962-445           |  |   |               |
| 337-16            | 50 kg | 800      |   | 337-160-200 |  |   | 962-446           |  |   |               |

#### Alternative to plastic container:

Wooden boxes for individual weights. For more details on this, please see page 206.

### Class F2 Test weights, stainless steel, stackable

Test weight material: finely turned stainless steel

Box material: Wood



| Test weight |       |          | +          | Container   |  | + | DAkkS certificate |  | = | Package price |
|-------------|-------|----------|------------|-------------|--|---|-------------------|--|---|---------------|
| KERN        |       | Tol ± mg | Dim. Ø x H | KERN        |  |   | KERN              |  |   | KERN          |
| 337-141     | 10 kg | 160      | 137x132 mm | 337-141-200 |  |   | 962-444           |  |   |               |
| 337-151     | 20 kg | 300      | 137x217 mm | 337-151-200 |  |   | 962-445           |  |   |               |
| 337-161     | 50 kg | 800      | 198x250 mm | 337-161-200 |  |   | 962-446           |  |   |               |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.