



PGR1200UK Large Pharmacy Refrigerator Glass Door 1200L

PGR1200UK

- ✓ 2°C to 8°C temperature range
- ✓ External digital controller and temperature display showing air and core temperatures
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ Temperature monitoring software records air and load temperatures constantly and also holds a full 7 day alarm history
- ✓ USB data download
- ✓ Power failure alarm

SPECIFICATION

Dimensions	1300mm W x 1995mm H x 840mm D
Weight	210kg
Product name	Large Pharmacy Refrigerator Glass Door 1200L PGR1200UK
GDPA code	PGR1200UK
SKU code	UK 444440366, EU 444442592
External dimensions	1300mm W x 1995mm H x 840mm D
Internal dimensions	1220mm W x 1386mm H x 653mm D
Packed dimensions	1480mm W x 2196mm H x 990mm D
Product weight	210kg
Packed weight	230kg
Gross capacity	1200L
External finish	Stainless steel

Colour	White/green
Internal finish	Stainless steel
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Glass
Door finish	Painted stainless steel
Door style	Hinged
Open door depth	1145mm
Replacement door seal	Y
Open door sensor	Y
Lock	Y
Adjustable feet	Fitted castors
Temperature range	2°C to 8°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	10
Shelf material	Wire
Number of shelf positions	40
Shelf size	645mm W x 523mm D
Interior light	Y
Auto or manual defrost	Auto
Fan	Y
Controller type	Digital, external
Climate class	ST
Refrigerant type	R600a
CFC free	Y
HFC free	Y
HCFC free	Y
Ammonia free	Y
UK 3 pin plug	Y
Voltage supply	230

Spark free	N
Auxiliary output	Y
Insulation thickness	60mm
Insulation type	Cyclopentane
Battery low alarm	Y (visual)
Air/load/defrost sensor fail alarm	Y (visual)
Compressor alarm	Y (visual)
Door open alarm	Y (visualand audio)
Air temperature alarm	Y (visualand audio)
Mains failure alarm	Y (visualand audio)
Antimicrobial handle	N
Optional drawer system available	N





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.