

# Control Plus

## Upright Control Plus Glass Door Refrigerator 353L

GDPA code **PPGR353UK**

SKU code **444410665**



### Key features

- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 5 year warranty as standard



#### New controller

New updated Controller system for improved user experience. Read, Record, Reset simply with the touch of one button, with no stress of accidentally resetting on the clear digital display.



#### Dual probe

Greater temperature accuracy given from two temperature probe readings. One probe measuring air temperature and one measuring load (or product) temperature.



#### Porthole

Factory fitted 16mm porthole enables quick and easy access for any secondary temperature monitoring probes.



#### Flexible internal storage

Multiple shelf locations along with a wire basket provides a flexible storage solution.



#### Anti-microbial handle

Powder coated to help prevent cross contamination by inhibiting the growth of microbes.

## Specification: Upright Control Plus Glass Door Refrigerator 353L

|                                     |                                             |
|-------------------------------------|---------------------------------------------|
| Gross Capacity (litres)             | 353L                                        |
| Net Capacity (litres)               | 349L                                        |
| Unit External Width (mm)            | 595mm                                       |
| Unit External Depth (mm)            | 640mm                                       |
| Unit External Height (mm)           | 1860mm                                      |
| Product Weight (kg)                 | 79kg                                        |
| Packed Width (mm)                   | 650mm                                       |
| Packed Depth (mm)                   | 727mm                                       |
| Packed Height (mm)                  | 1885mm                                      |
| Packed Weight (kg)                  | 85kg                                        |
| Country of Origin                   | UK                                          |
| Internal Dimensions (mm) (WxHxD)    | 480mm W x 1682mm H x 445mm D                |
| Temperature Range                   | +2°C to +8°C (set point tolerance + or - 2) |
| Colour                              | White/Green/Glass                           |
| External Finish                     | Painted aluminium                           |
| Internal Finish                     | HIPS*                                       |
| Type of plastic                     | SRL800 (HIPS/ABS mix)                       |
| Door Finish                         | Painted aluminium/glass                     |
| Door Type                           | Painted aluminium/glass                     |
| Door Style                          | Hinged right hand only, non reversible      |
| Door Liner                          | Flat                                        |
| Door Open Depth (mm)                | 1800mm                                      |
| Replacement Door Seal               | Yes                                         |
| Open Door Sensor                    | Yes                                         |
| Lock                                | Yes, 2 Keys                                 |
| Adjustable Feet (number & position) | 2 adjustable front feet, 2 rear rollers     |
| Porthole                            | Yes                                         |
| Porthole Size                       | 16 mm                                       |
| Porthole Location                   | Rear Bottom                                 |
| Temperature Display on Panel        | Yes                                         |
| Temperature Alarm                   | Yes                                         |
| Remote Alarm Contacts               | Yes                                         |
| Audible Alarms                      | Yes                                         |
| Visual Alarms                       | Yes                                         |
| Number of Shelves                   | 6                                           |
| Number of Shelf Positions           | 14                                          |
| Shelf Size (mm) (WxDxH)             | 385mm W x 6mm H x 500mm D                   |
| Shelf Material                      | Wire                                        |
| Number of Baskets                   | 1                                           |
| Basket Material                     | Wire                                        |
| Basket Size (mm) (WxDxH)            |                                             |
| Interior Light                      | Yes                                         |
| Light Type                          | N/A                                         |
| Light Activation                    | N/A                                         |
| Fan                                 | Yes                                         |
| Spark Free                          | No                                          |
| Defrost setting                     | Auto                                        |
| Controller Type                     | Digital, External                           |
| Number of Probes                    | 2                                           |
| Probe Type                          | Air & Load (Glycol Solution)                |
| Ambient Temperature Range           | 10 to 32                                    |
| Refrigerant Type                    | R600a                                       |
| CFC Free                            | Yes                                         |
| HFC Free                            | Yes                                         |
| HCFC Free                           | Yes                                         |
| Ammonia Free                        | Yes                                         |
| UK 3 Pin Plug                       | Yes                                         |
| Other Plug Types Available          | Yes                                         |
| Voltage Supply                      | 220-240                                     |
| Insulation Thickness (mm)           | 60mm                                        |
| Insulation Type                     | High Density Polyurethane Foam              |
| Blowing Agent                       | Cyclopentane                                |
| Anti Microbial Handle               | Yes                                         |
| Wall Brackets                       | No                                          |
| Optional Drawer System              | No                                          |
| Optional Digital Lock               | No                                          |
| Other Door Option                   | Solid PPSR353UK - 444410666                 |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.