

# Pharmacy Plus Bluetooth Enabled Refrigerator

## 120L Pharmacy Plus Undercounter Fridge with Solid Door

Product code PPSR120BT

SKU code 444411883




### Key features

- ✓ Temperature range of +2°C to +8°C
- ✓ Audible and visual temperature alarms
- ✓ Solid and glass door options
- ✓ Integrated lock, supplied with 2 keys for added security



#### Bluetooth Connected Controller

Our controller is Bluetooth enabled and allows for simple, contactless monitoring of temperature and alarm history, via the official Lec Medical app. Battery back-up, provides continuous temperature monitoring in the event of power outage.



#### Lec Medical Smartphone App

Contactless local monitoring of air and load temperatures via smart phone app. Downloadable temperature report and alarm history.



#### Intelligent Fan Management

Smart fan optimisation aides in temperature stability following door openings, minimising fan activity when door is open to reduce the rate of temperature increase.



#### Dual Probe Temperature Monitoring

Our dual probe system offers greater temperature accuracy within your appliance, simultaneously measuring both air and product temperatures.



#### Flexible Internal Storage

Multiple shelf locations provides flexible storage solutions.

Specification: 120L Pharmacy Plus Undercounter Fridge with Solid Door

| Product Information          |                       | Door Specification                  |                                |
|------------------------------|-----------------------|-------------------------------------|--------------------------------|
| Brand                        | Lec Medical           | Door Finish                         | Painted steel                  |
| Category                     | Pharmacy              | Door Type                           | Solid                          |
| Sub-Brand                    | Pharmacy Plus         | Door Style                          | Hinged reversible              |
| SKU Code                     | 444411883             | Door Liner                          | Flat                           |
| EAN Code                     | 5052263118836         | Number of Door Trays                | n/a                            |
| GDMA Code                    | PPSRI20BT             | Number of Door Tray Positions       | n/a                            |
| Description                  | Under Counter         | Door Open Depth (mm)                | 1035mm                         |
|                              |                       | Replacement Door Seal               | Yes                            |
| Product Features             |                       | Technical Information               |                                |
| Gross Capacity (litres)      | 120L                  | Noise Level (dB)                    | 43.0dB                         |
| Net Capacity (litres)        | 112L                  | Colour                              | White/green                    |
| Temperature Range            | +2°C to +8°C          | External Finish                     | Painted steel                  |
| Open Door Sensor             | Yes                   | Internal Finish                     | HIPS                           |
| Key Lock                     | Yes, 2 Keys           | Type of plastic                     | SRL800 (HIPS/ABS mix)          |
| Digital Lock                 | No                    | Ambient Temperature Range           | 16°C to 32°C                   |
| Porthole                     | Yes                   | Refrigerant Type / Weight (g)       | R600a/26g                      |
| Porthole Size (Ø Diameter)   | 16mm                  | CFC Free                            | Yes                            |
| Porthole Location            | Rear Bottom (Right)   | HFC Free                            | Yes                            |
| Controller Type              | Digital BT, External  | HCFC Free                           | Yes                            |
| Temperature Display on Panel | Yes                   | Ammonia Free                        | Yes                            |
| Temperature Alarm            | Yes                   | Insulation Thickness (mm)           | 40mm                           |
| Remote Alarm Contacts        | Yes                   | Insulation Type                     | High Density Polyurethane Foam |
| Audible Alarms               | Yes                   | Blowing Agent                       | Cyclopentane                   |
| Visual Alarms                | Yes                   | Certification                       | CE/UKCA                        |
| Interior Light               | Yes                   | ATEX Zone Classification            | n/a                            |
| Light Type                   | LED                   |                                     |                                |
| Light Activation             | Door Switch / App     |                                     |                                |
| Fan - Type                   | Yes - Axial           |                                     |                                |
| Spark Free                   | No                    |                                     |                                |
| Defrost setting              | Automatic             |                                     |                                |
| Battery Back Up/Hours        | Yes/48hrs             |                                     |                                |
| Bluetooth                    | Yes                   |                                     |                                |
| NFC                          | Yes                   |                                     |                                |
| WIFI                         | No                    |                                     |                                |
| GSM-4G                       | No                    |                                     |                                |
| Number of Probes             | 2                     |                                     |                                |
| Probe Type                   | Air & Load            |                                     |                                |
| Anti Microbial Handle        | Yes                   |                                     |                                |
| Internal Specification       |                       | Power/Installation                  |                                |
| Number of Shelves            | 3                     | UK 3 Pin Plug Power Cord            | Yes                            |
| Number of Shelf Positions    | 7                     | EU Schuko 2 pin plug Power Cord     | Yes                            |
| Shelf Size 1 (mm) (WxDxH)    | 400mm x 390mm x 6mm   | AU Type I Plug Power Cord           | No                             |
| Shelf Size 2 (mm) (WxDxH)    | 400mm x 285mm x 6mm   | Voltage Supply (V)                  | 220-240V                       |
| Shelf Material               | Polythene Coated Wire | Frequency (Hz)                      | 50/60Hz                        |
| Number of Drawers            | n/a                   | Max Power (W)                       | 86.0W                          |
| Number of Drawer Positions   | n/a                   | Max Current Rating (A)              | 0.60A                          |
| Drawer Size 1 (mm) (WxDxH)   | n/a                   | Energy Consumption (kWh/24h)        | 0.37kWh/24h                    |
| Drawer Size 2 (mm) (WxDxH)   | n/a                   | Energy Consumption AE (kWh)         | 135kWh                         |
| Drawer Material              | n/a                   | Heat Output Estimated at 22°C (kJ)  | 985kJ/h                        |
| Number of Baskets            | 1                     | Adjustable Feet (Number & position) | 2 adjustable (front)           |
| Basket Material              | Polythene Coated Wire |                                     |                                |
| Basket Size (mm) (WxDxH)     | 375mm x 295mm x195mm  |                                     |                                |
|                              |                       | Weights & Dimensions                |                                |
|                              |                       | Unit External Width (mm)            | 475mm                          |
|                              |                       | Unit External Depth (mm)            | 595mm                          |
|                              |                       | Unit External Height (mm)           | 845mm                          |
|                              |                       | Unit External Depth Inc Handle (mm) | 635mm                          |
|                              |                       | Internal Dimensions (mm) (WxDxH)    | 385mm x 495mm x 700mm          |
|                              |                       | Product Weight (kg)                 | 30.14Kg                        |
|                              |                       | Packed Width (mm)                   | 505mm                          |
|                              |                       | Packed Depth (mm)                   | 660mm                          |
|                              |                       | Packed Height (mm)                  | 865mm                          |
|                              |                       | Packed Weight (kg)                  | 33.96Kg                        |
|                              |                       | Country of Origin                   | China                          |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.