

Bluetooth Enabled Pharmacy Plus Range

47L Bluetooth Pharmacy Countertop Fridge Glass Door

GDPA code PPGR47BT-UK

SKU code 444411306




Key features

- ✓ Local Monitoring and Alarm History via Bluetooth
- ✓ Power Failure Alarm with battery back up. Rechargeable battery
- ✓ Solid and glass door options
- ✓ Intelligent Fan Management for Enhanced Performance



New Controller

Offers Local temperature monitoring and alarm history via Bluetooth. Power failure alarm with battery backup, rechargeable battery and intelligent fan management for enhanced performance



New Bluetooth Enabled App

Accessibility features include contactless local monitoring of air and load temperatures. Downloadable temperature history and alarm history.



Flexible internal storage

Multiple shelf locations along with a wire basket provides a flexible storage solution.



Anti-microbial handle

Powder coated to help prevent cross contamination by inhibiting the growth of microbes.



Dual probe

Greater temperature accuracy given from two temperature probe readings. One probe measuring air temperature and one measuring load (or product) temperature.

Specification: Upright 41L Bluetooth Pharmacy Countertop Fridge Glass Door

Gross Capacity (litres)	47L
Net Capacity (litres)	30L
Unit External Width (mm)	475mm
Unit External Depth (mm)	445mm
Unit External Height (mm)	565mm
Product Weight (kg)	17.5kg
Packed Width (mm)	490mm
Packed Depth (mm)	530mm
Packed Height (mm)	620mm
Packed Weight (kg)	21.5kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	403mm W x 354mm D x 425mm H
Temperature Range	2°C to 8°C
Colour	White/green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	PVC/Glass
Door Type	Glass
Door Style	Hinged, reversible
Door Liner	Flat
Door Open Depth (mm)	905mm
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	2 adjustable (front)
Porthole	Yes
Porthole Size	16 mm
Porthole Location	Rear LH Corner
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	1
Number of Shelf Positions	3
Shelf Size (mm) (WxDxH)	416mm W x 200mm D x 6mm H
Shelf Material	Polythene Coated Wire
Number of Baskets	1
Basket Material	Polythene Coated Wire
Basket Size (mm) (WxDxH)	396mm W x 165mm D x 150mm H
Interior Light	Yes
Light Type	LED
Light Activation	Door Swith / Controller
Fan	Yes - Axial
Spark Free	No
Defrost setting	Auto (Off Cycle)
Controller Type	Digital BT, external
Number of Probes	2
Probe Type	Air & Load (Glycol Solution)
Ambient Temperature Range	10°C to 32°C
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	220 -240V
Insulation Thickness (mm)	40mm
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Wall Brackets	Yes
Optional Drawer System	No
Optional Digital Lock	No
Other Door Option	Solid PPSR47UK - 444410570
Battery Backup	Yes



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.